

# BLICKMAN



## ***Exam Stools***

*Providing ultimate comfort and function*

**Includes:**

- Spindle Adjusted Stools
- Pneumatic Stools
- Parts & Accessories



# THE BLICKMAN DIFFERENCE

## EXAM STOOLS

- Our durable, high quality chrome and stainless steel exam stools come with 3" thick, dense foam 15" diameter cushion for ultimate comfort.
- Our upholstery is washable vinyl and is available in an array of 75 fashion colors that meet specific design requirements - see upholstery catalog for details.
- For applications that require the longevity that stainless steel provides, pick from our revolving stools, or for a more modern look, choose from four or five leg models with spindle adjusted or pneumatic lift mechanisms.
- For back support, choose one of our models with a back rest.
- Our three styles of locking casters allow for easy mobility and safety.

## Key Product Features



# Exam Stools

Spindle

## 4 LEG SPINDLE EXAM STOOLS

**Description:** We offer two four leg spindle adjusted exam stools with the option of 2" casters or rubber crutch tips. Choose our Passaic stool and sit comfortably on a 15" diameter, 3" thick polyfoam padded seat upholstered in vinyl or choose our Clifton stool for an easy-to-clean stainless steel seat with a non-slip vinyl insert. Add a removable pad slipcover to your Clifton stool for a more comfortable option.



### Spindle

Adjusts 12"  
**Rubber Tips:** 19" - 31"  
**Casters:** 21" - 33"



### Name

### CLIFTON

### PASSAIC

#### Standard

#### With Pad

#### Standard

#### Casters

Model #	7745SS-C	7745SS-PC	7714-PSSC
Item #	1027745001	1027745003	1027714125
Shipping Weight	17 lbs.	17 lbs.	17 lbs.

#### Crutch Tips

Model #	7745SS	7745SS-P	7714-PSS
Item #	1027745000	1027745002	1027714025
Shipping Weight	16 lbs.	16 lbs.	16 lbs.



**Options:** Available in standard black. Consult factory for requests in other colors for large order quantities

# Exam Stools

Spindle

## 5 LEG SPINDLE EXAM STOOL

**Description:** Our five leg spindle exam stool rests on a 24" diameter base. Choose between a black composite base or an aluminum base with Wall Saver™ tips on the ends to protect against scuff marks and offer a textured surface for foot grip. Sit comfortably on its 15" diameter, vinyl upholstered seat with 3" thick polyfoam padding. This model is available with or without a back.



**Spindle**

Adjusts 5 7/8"  
**Height:** 20 7/8" - 26 3/4"



Name	5 Leg Spindle Stool	
	Aluminum Base	5 Leg Spindle Stool with backrest
Model #	1113	1113WB
Item #	1021113000	1021113200
Shipping Weight	22 lbs.	26 lbs.
	Black Composite Base	
Model #	1113BLK	1113WB-BLK
Item #	1021113100	1021113300
Shipping Weight	22 lbs.	26 lbs.

**Options:** Foot ring, 75 different vinyl colors - see upholstery catalog

# Exam Stools

Pneumatic

## FOOT OPERATED EXAM STOOL

**Description:** Our five-leg foot operated exam stool have a push ring under a 24" diameter aluminum base with Wall Saver™ tips on the ends to protect against scuff marks and offer a textured surface for foot grip. Sit comfortably on its 15" diameter, vinyl upholstered seat with 3" thick polyfoam padding. This model is available with or without a back.



**Foot Pedal**

Adjusts 7"  
Height: 18 3/8" - 25 3/8"



Name	Foot Operated Stool	Foot Operated Stool with Backrest
	<b>Aluminum Base</b>	
Model #	1205	1206
Item #	1041205000	1041206000
Shipping Weight	22 lbs.	28 lbs.

**Options:** Foot ring, 75 different vinyl colors - see upholstery catalog

# Exam Stools

Pneumatic

## HAND OPERATED EXAM STOOLS

**Description:** Our five-leg hand operated exam stool rests on a 24" diameter base. Choose between a black composite base or an aluminum base with Wall Saver™ tips on the ends to protect against scuff marks and offer a textured surface for foot grip. Sit comfortably on its 15" diameter, vinyl upholstered seat with 3" thick polyfoam padding. This model is available with or without a back.

### Standard Stools



#### Hand Lift

Adjusts 7"  
**Height:** 18 3/8" - 25 3/8"



Name	Hand Operated Stool	Hand Operated Stool with Backrest
	<b>Aluminum Base</b>	
Model #	1210	1210WB
Item #	1041210000	1041210200
Shipping Weight	22 lbs.	26 lbs.
	<b>Black Composite Base</b>	
Model #	1210BLK	1210WB-BLK
Item #	1041210100	1041210300
Shipping Weight	22 lbs.	26 lbs.

**Options:** Foot ring, 75 different vinyl colors - see upholstery catalog



# Exam Stools

Pneumatic

## HAND OPERATED EXAM STOOLS

### Economy Stools

**Description:** Our economy hand operated exam stool rests on a 24" diameter black composite base with black plastic twin wheel casters. Sit in comfort on a 15" diameter, vinyl upholstered seat with 3" thick polyfoam padding. this model is available with and without a back.

Adjusts 7" from 18 3/8" to 25 3/8" off the ground



Name	Economy Stool	
	Without Back	With Back
Model #	1200	1201
Item #	1041200100	1041201300
Shipping Weight	22 lbs.	26 lbs.

**Options:** Foot ring, 73 different vinyl colors - see upholstery catalog

### Deluxe Stools

**Description:** Our deluxe hand operated exam stool comes standard with a foot ring for leg-joint comfort. These stools rest on a 24" diameter black composite or aluminum base, which comes with Wall Saver™ tips to protect against scuff marks and offer a textured surface for foot grip. Both bases travel easily on 2" hard rubber casters. Sit in comfort on a 15" diameter, vinyl upholstered seat with 3" thick polyfoam padding with a backrest.

Adjusts 7" from 21 3/8" to 28 3/8" off the ground



Name	Deluxe Stool	
	Aluminum Base	Black Composite Base
Model #	1212	1212BLK
Item #	1041212000	1041212100
Shipping Weight	29 lbs.	29 lbs.

**Options:** 75 different vinyl colors - see upholstery catalog

**Foot Ring as an Accessory:** Compatible on all 5 leg stools

Part #	Description	Shipping Weight
7901243000	Foot Ring for Leg Joint Comfort	6 lbs.





# Exam Stools

## Construction Specifications

### Construction Specifications

Standard Dimensions	<b>Seats:</b> 15" diameter <b>Polyfoam padding:</b> 3" thick Height Range: 18" - 33"
Materials	<b>Bases:</b> Aluminum or black composite <b>Seats:</b> 75 colors of upholstered vinyl covering polyfoam padding. Clifton/Passaic Stool: 304 stainless steel
Gauge	Clifton Stool Legs: 18GA; Clifton Stool Seat: 14 GA Passaic Stool Legs : 18 GA
Welding	Clifton/Passaic Stool: All welding by heliarc method and exposed welds are ground smooth
Finish	Clifton/Passaic Stool: Exposed surfaces are free from pit marks, weld seams, scale and are polished with #4 polish to a satin finish
Casters	2" hard rubber casters or 2" twin wheel plastic casters
Notes	<ul style="list-style-type: none"> <li>• MR conditional stool available. See MR conditional catalog</li> <li>• For all vinyl colors, please see our upholstery catalog</li> </ul>

### Shipping Information

All shipping weights are estimates and may be changed without prior notice. Specifications may be improved and dimensions changed without prior notice. Please contact customer service for estimated shipping costs at [info@blickman.com](mailto:info@blickman.com).

Toll Free: (800) 247-5070  
Local: (973) 330-0557



[www.blickman.com](http://www.blickman.com)  
[info@blickman.com](mailto:info@blickman.com)

1099 Wall Street West, Suite 250A, Lyndhurst, NJ 07071