



EXAM/PATIENT ROOM EQUIPMENT

Engineered for a Lifetime — Our Reputation is Stainless



CHROME EXAM STOOLS
FOOT STOOLS
WASTE RECEPTACLES
MAYO STANDS

UTILITY CARTS
PRIVACY SCREENS
OVERBED TABLES



CHROME



1111

MODEL 1111 CHROME FOUR-LEG ADJUSTABLE EXAM STOOL is designed for smooth operation and comfort with a chrome shaft and nylon bushing that rides quietly and cleanly during adjustment. Features a 15"-diameter seat, 3" thick polyfoam padding on durable vinyl upholstery. Stool travels on 2", single-wheel, hard-rubber casters. Choose from 12 upholstery colors.

Model No.	Feature	Part No.	Height Adjustment	Ship Wt. (lb) est.
1111	round foot ring	1021111000	17-1/2"-27-1/2"	15

MODEL 1113BLK is our five-leg adjustable exam stool that rests on a black composite base for a contemporary look. Sit in comfort with a 15"-diameter seat, 3" thick polyfoam padding on durable vinyl upholstery. Stool travels on 2", single-wheel, hard-rubber casters. Choose from 12 upholstery colors.

Model No.	Feature	Part No.	Height Adjustment	Ship Wt. (lb) est.
1113BLK	black base	1021113100	20-1/2"-26-1/4"	22
1113WB/BLK	backrest black base	1021113300	20-1/2"-26-1/4"	25



1113BLK



1210WB/BLK

MODEL 1210WB/BLK FIVE-LEG PNEUMATIC STOOL offers smooth height adjustment with the touch of the lever on the Supra Pneumatic Lift™. Sit comfortably on the wide 15"-diameter vinyl seat cushion supported by 3" of dense foam padding. The black composite base provides a contemporary look. Two-inch, single-wheel, hard-rubber casters roll easily on carpet or floor. Choose from 12 upholstery colors.

Model No.	Feature	Part No.	Height Adjustment	Ship Wt. (lb) est.
1210BLK	black base backrest	1041210100	22"-29"	25
1210WB/BLK	black base backrest	1041210300	22"-29"	25

CHROME



MODELS 1250, 1251 FOOT STOOLS are constructed of chrome and have 17 weld points to enhance strength and durability. These stools feature a non-slip rubber mat and rubber tips for added comfort. Add the optional handrail – Model 1251.

Model No.	Part No.	Dimensions	Ship Wt. (lb) est.
1250	1011250000	11-3/8"D 14-1/4"W 9"H	8
1251(handrail)	1011251000	11-3/8"D 14-1/4"W 9"H, (35-3/4"HR)	10

Option Shown: Handrail attachment (model 1251).



MODELS 2820 AND 2830 WASTE RECEPTACLES feature a white painted step-on can with a non-slip pedal and rubber lid guard for secure closure. Black poly pail insert.

Model No.	Feature	Part No.	Dimensions	Ship Wt. (lb) est.
2820	20-Qt.	0832820000	11-1/2"D 11"W 17"H	12
2830	32-Qt.	0832830000	12-1/4"D 11-1/4"W 18-1/4"H	15



STAINLESS



2840SS

MODEL 2840SS STAINLESS STEEL WASTE RECEPTACLE provides a stainless steel step-on can with a quiet, recessed pedal and rubber lid guard for secure closure. Includes a removable poly insert.

Model No.	Feature	Part No.	Dimensions	Ship Wt. (lb) est.
2840SS	20-Qt	0832840000	11-3/4"D 9-3/4"W 17-3/4"H	10

CHROME

MODEL 1501 INSTRUMENT STAND features a five-leg cast aluminum base with wall savers. A Tru- Loc™ friction knob secures the height adjustment of the stainless steel tray. Travels easily on 2", twin-wheel, hooded casters.

Model No.	Feature	Part No.	Height Adjustment	Tray Size	Ship Wt. (lb) est.
1501	Five-leg base	0661501000	29" to 49"	12-5/8"x19-1/8"	12




1501



1510

DOUBLE-POST AND SINGLE-POST CHROME MAYO STANDS adjust easily with the Tru- Loc™ friction knob to ensure secure positioning of the stainless steel tray. Travels easily on 2" twin wheel casters (1510) or hard rubber casters (1515).

Model No.	Feature	Part No.	Height Adjustment	Tray Size	Ship Wt. (lb) est.
1510	Double-post	0661510000	34-1/2" to 54"	12-5/8" x 19-1/8"	19
1515	Single-post	0661515000	31-1/2" to 51"	12-5/8" x 19-1/8"	12





CHROME



2440

MODELS 2440, 2441, 2442 FOLDING UTILITY CARTS are versatile enough for any application such as movement of mail, patient records, and laboratory materials. Choose any combination of 6" or 12" deep poly-coated wire baskets to complete your cart. The metal frame is designed for easy assembly and folds for storage. Rubber stoppers support and protect the poly-coated wire racks and provide a quiet ride on four, three-inch hard rubber casters.

Model No.	Feature	Part No.	Dimensions	Ship Wt. (lb) est.
2440	folding cart	2422440000	19-1/2"D 29-3/4"W (floor space) 35" H	21
2441	6" wire basket	2422441000	14-1/4"D 24-1/2"W 6"H	5
2442	12" wire basket	2422442000	14-1/2"D 24-3/4"W 12"H	6

MODEL 2900 AND 2904 PATIENT PRIVACY SCREENS are fabricated from sturdy anodized aluminum with stainless steel hinges. The panels are self-extinguishing white vinyl and fold easily for compact storage. All models meet flammability requirements of California Bulletin 117.



2900



2904

Model No.	Feature	Part No.	Dimensions	Ship Wt. (lb) est.
2900	Three-Panel	2822900000	51"L (extended) x 19"W (folded) x 67-1/2"H	9
2904	Three-Panel with casters	2822904000	78-1/4"L (extended) x 26-1/2" W (folded) x 70-1/2" H	19



3400

MODEL 3400 OVERBED TABLE moves up or down with a simple one-touch adjustment. Choose an attractive mar-resistant laminate top in wood grain. Security features include a built-in column lock and four-point base. The table travels smoothly on four, two-inch hooded casters.

Model No.	Feature	Part No.	Dimensions	Ship Wt. (lb) est.
3400	wood grain	2833400000	Base-15-1/2"D x 26-1/2"W Top-15"x30"H	27


EXAM ROOM EQUIPMENT

Model No.	Part No.	Description	Ship Wt. (lb) est.
1111	1021111000	4-Leg Adjustable Stool with Round Foot Ring	15*
1113BLK	1021113100	5-Leg Adjustable Stool with Black Base	22*
1113WB/BLK	1021113300	5-Leg Adjustable Stool with Black Base & Backrest	25*
1210BLK	1041210100	5-Leg Pneumatic Stool with Backrest	25*
1210WB/BLK	1041210300	5-Leg Pneumatic Stool with Black Base & Backrest	25*
1250	1011250000	Foot Stool	8*
1251	1011251000	1250 with Handrail. Overall ht. 35-3/4"	10*
2820	0832820000	Step-On Can 20-Qt.	12*
2830	0832830000	Step-On Can 32-Qt.	15*
2840SS	0832840000	Stainless Step-On Can 20-Qt.	10*
1501	0661501000	5-Leg Instrument Stand with s/s tray	12*
1510	0661510000	Double-Post Mayo Stand with s/s tray	19*
1515	0661515000	Single-Post Mayo Stand (U-base)	12*
2440	2422440000	Folding Utility Cart	21*
2441	2422441000	6" Wire Basket	5*
2442	2422442000	12" Wire Basket	6*

PATIENT ROOM EQUIPMENT

Model No.	Part No.	Description	Ship Wt. (lb) est.
2900	2822900000	Patient Privacy Screens- 3- Panel	9*
2904	2822904000	Patient Privacy Screens- 3 Panel with Casters	19*
3400	2833400000	Overbed Table, Wood Grain	27*

* Items can be shipped via FedEx® Ground delivery service.



39 Robinson Road Lodi, NJ 07644 • 800-247-5070 • Fax: 201-909-0832
e-mail: info@blickman.com • www.blickman.com

All shipping weights are estimates only and may be changed without prior notice. Specifications may be improved and dimensions changed without prior notice. Unless otherwise noted, all products displayed or listed are deemed, MR-'Unsafe,' in accordance with ASTM Standard F2052-06.

FedEx service marks used by permission.

Copyright © 2006 Blickman, Inc.