



Note: Photos depict the bassinets with **optional** infant basket and pad, however, please order them separately to avoid disappointment.



MODEL 8044SS BASSINET has two shelves welded to the inverted “U” frame that measure 8” and 23” from the floor. Upper shelf has a four-sided guardrail to protect supplies. All welded construction. Four-inch, ball-bearing, swivel casters, two with brakes, transport baby comfortably and safely.

Model No.	Part No.	Feature	Dimensions (overall)	Ship Wt. (lb) est.
8044SS	1118044000	two shelves	31"W 17-3/4"D 38-1/2"H	100

Option: H-brace in lieu of lower shelf.



MODEL 8047SS BASSINET provides a storage drawer and extra shelf 8” from the floor. The drawer is mounted to separate channels instead of the basket rods. Drawer dimensions: 21-1/2”W 16”D 7-1/2”H. All welded construction. Four-inch, ball-bearing, swivel casters, two with brakes, transport baby comfortably and safely. Counter-weighted drawer aids stability.

Model No.	Part No.	Feature	Dimensions (overall)	Ship Wt. (lb) est.
8047SS	1118047000	drawer & shelf	31"W 17-3/4"D 38-1/2"H	100



MODEL 8048SS BASSINET provides a two-door cabinet welded to the inverted “U” frame design. The cabinet has two 17” deep storage compartments; the right side has a removable shelf. Cabinet dimensions: 29”W 16-3/4”D 15”H. All welded construction. Four-inch, ball-bearing, swivel casters, two with brakes, transport baby comfortably and safely.

Model No.	Part No.	Feature	Dimensions (overall)	Ship Wt. (lb) est.
8048SS	1118048000	2-door cabinet	31"W 17-3/4"D 38-1/2"H	100



Note: Photos depict the bassinets with **optional** infant basket and pad, however, please order them separately to avoid disappointment.



MODEL 8049SS BASSINET offers an open cabinet with two storage compartments, one with shelf. Cabinet dimensions: 29"W 17"D 15"H. All welded construction. Four-inch, ball-bearing, swivel casters, two with brakes, transport baby comfortably and safely.

Model No.	Part No.	Feature	Dimensions (overall)	Ship Wt. (lb) est.
8049SS	1118049000	Open cabinet	31"W 17-3/4"D 38-1/2"H	100


BASSINET ACCESSORIES

These items are not included with bassinets; please order separately to avoid disappointment.



MODEL 3070P INFANT BASKET - This seamless basket of clear molded plastic has a cardholder for easy identification. Blickman bassinets utilize specially formed support rods so basket can be set in the flat or trendelenburg position.

Model No.	Part No.	Dimensions (overall)	Ship Wt. (lb) est.
3070P	9841307001	31"L 16"W 9-1/4"H (head end)	10



MODEL 3069MP MATTRESS PAD - Two-inch thickness provides comfort to newborn. Staph check cover over foam padding. Flame retardant.

Model No.	Part No.	Dimensions (overall)	Ship Wt. (lb) est.
3069MP	8449306901	12-1/2"L 26-1/2"W 2"H	2



PRODUCT LISTING

BASSINETS

Model No.	Part No.	Description	Dimensions (Overall)	Ship Wt. (lb) est.
8044SS	1118044000	Bassinet with two shelves	31"W (78.74) x 17-3/4"D (45.09) x 38-1/2"H (97.79)	100
8047SS	1118047000	Bassinet with storage drawer and shelf	31"W (78.74) x 17-3/4"D (45.09) x 38-1/2"H (97.79)	100
8048SS	1118048000	Bassinet with two-door cabinet	31"W (78.74) x 17-3/4"D (45.09) x 38-1/2"H (97.79)	100
8049SS	1118049000	Bassinet with open cabinet	31"W (78.74) x 17-3/4"D (45.09) x 38-1/2"H (97.79)	100

BASSINET ACCESSORIES

INFANT BASKET AND PAD

3070P	9841307001	Infant Basket	31"L (78.74) x 16"W (40.64) x 9-1/4"H (23.50)	10*
3069MP	8449306901	Mattress Pad	12-1/2"L (31.75) x 26-1/2"W (67.31) x 2"H (5.08)	2*

* Items can be shipped via FedEx® Ground delivery service.



39 Robinson Road Lodi, NJ 07644 • 800-247-5070 • Fax: 201-909-0832
e-mail: info@blickman.com • www.blickman.com

All shipping weights are estimates only and may be changed without prior notice. Specifications may be improved and dimensions changed without prior notice. Unless otherwise noted, all products displayed or listed are deemed, MR-'Unsafe,' in accordance with ASTM Standard F2052-06.

FedEx service marks used by permission.

Copyright © 2006 Blickman, Inc.