

**BLICKMAN
BUILT™**

HAMPERS/RECEPTACLES

Engineered for a Lifetime — Our Reputation is Stainless



Part No.	Description	Qty./Box
1. 760	Body	1 PA
2. 760	Top	1 PA
3. 760	Bracket Assy	1 PA
4. 760	Bracket Cap	1 PA
5. 760	Roller Assembly	1 PA
6. 870000000	Capset Pk Assy	2 PA
7. 870000000	Capset Pk Assy	2 PA

BUCKMAN HEALTH INDUSTRIES, INC.
 1000 N.J.
 N.A.
 100 EAST 100
 PHILADELPHIA, PA 19106
 TEL: 215-381-1000
 FAX: 215-381-1001
 WWW.BUCKMANHEALTH.COM
 © 2000 BUCKMAN HEALTH INDUSTRIES, INC.



HAMPERS
RECEPTACLES



CHROME



2010

Keep linens under cover with the foot-operated **MODEL 2010 STEP-ON HAMPER**. Just press the non-slip pedal to raise the white poly-coated steel lid; locking device keeps it in open position as needed. Supports 36-42 gallon multi-purpose bags with a ring to keep bags in place. Includes hamper bag storage pouch (2012). Glides effortlessly on 3" hard wheel casters, two with locks.

Model No.	Part No.	Dimensions	Ship Wt. (lb) est.
2010	0962010000	LID: 18-5/8"W 19-1/2"D 37-3/4"H	25

Option: Hamper bags



2018

MODEL 2018 and **2025 ROUND HAMPERS** offer an 18" or 25" diameter ring constructed from 1" chrome steel tubing. Solid steel bag rest supports hamper bags. Travels easily on 2" hard rubber casters.

Model No.	Part No.	Dimensions	Ship Wt. (lb) est.
2018	0962018000	18-5/8"D 35"H	11
2025	0962025000	25"D 35"H	12

Option: Hamper bags



2118

MODEL 2118 TILT-TOP FOOT-OPERATED HAMPER is spring-loaded for easy opening and closing and offers an attractive stainless steel lid. Non-slip pedal allows for hands-free operation. Solid steel bag rest supports hamper bags. The hamper travels easily on four, 2" hard rubber casters.

Model No.	Part No.	Dimensions	Ship Wt. (lb) est.
2118	0962118000	19-5/8"D 36"H	21

Option: Hamper bags





2230

MODEL 2230 FOLDING HAMPER saves space and adjusts to accommodate three hamper bag sizes: 18", 22" or 25". Steel frame bag rest supports hamper bags. Simple to operate. This unit travels easily on four, 2" hard rubber casters.

Model No.	Part No.	Dimensions		Ship Wt. (lb) est.
2230	0962230000	36-1/2"H	19"W	18
		40-1/2"H	19"W	
		42-1/4"H	19"W	

Option: Hamper bags

STAINLESS



2301SS

MODEL 2301SS TRIANGULAR SPACE-SAVER HAMPER offers hands-free loading and unloading with a convenient non-slip pedal. All stainless steel construction. Bumper guards on stainless steel lid and legs protect equipment and walls. Solid shelf on bottom provides extra security with heavy loads. This hamper travels effortlessly on three, 3" hard rubber swivel casters.

Model No.	Part No.	Dimensions		Ship Wt. (lb) est.
2301SS	0962301000	18"D	32"H	25

Option: Hamper bags

STAINLESS



MODELS 7773 HAMPERS have round, tubular rings that are seamlessly welded to three legs. A solid, round platform at the bottom acts as a bracing and protects the hamper bag from the floor. For Model 7773SS, an optional, reinforced, dome-shaped lid of high-density, white polyethylene is recommended to control airborne pathogens. MODEL 7773SS-LP comes standard with a plastic lid, while MODEL 7773SS-LPF is also foot-operated.

MODELS 7774SS ROUND HAMPER offers an 18-5/8" diameter ring, constructed from 1" stainless-steel tubing. Solid steel bag-rest supports hamper bags. MODEL 7774SS-LF TILT-TOP FOOT-OPERATED HAMPER is spring-loaded for easy opening and closing, and offers an attractive, 20"-diameter, stainless-steel lid. Non-slip pedal allows for hands-free operation. Solid steel bag rest supports hamper bags. Both hampers travel easily on 2" hard rubber casters.

Model No.	Part No.	Lid	Dimensions (overall)		Ship Wt. (lb) est.
7773SS	0927773000	none	26"D	35"H	20
7773SS-LP	0927773101	plastic	28"D	36.5"H	35
7773SS-LPF	0927773201	plastic; foot-op	28"D	36.5"H	55
7774SS	0927774000	none	18-5/8"D	35"H	11
7774SS-LF	0927774201	stainless-steel; foot-op	19-5/8"D	36"H	21

OPTIONAL HAMPER BAGS:

Model No.	Part No.	Description
HB18	9839777201	18" Diameter Bags
HB44	9839775202	44" Diameter Bags

OPTIONAL PLASTIC LID KIT:

Model No.	Part No.	Description
PLK-26	0927773548	Plastic Lid Kit-26"



STAINLESS



7766SS

MODEL 7766SS (LENOX) KICKBUCKET has a continuous encircling rubber bumper with stainless steel cross braces that provide excellent protection for receptacle. Constructed from heavy gauge stainless steel, this 13 qt. pail fits perfectly into a tip-resistant base protected by rubber tips for smooth handling and noise reduction. Unit travels effortlessly on 2" swivel casters.

Model No.	Part No.	Dimensions	Capacity	Ship Wt. (lb) est.
7766SS	0817766000	16"D 12-1/2"H	13 qt.	12



8766SS

MODEL 8766SS SPONGE RECEPTACLE offers the same high quality construction as the Lenox Kickbucket but accommodates a round basin with an 8-1/2 quart capacity. Unit travels effortlessly on 2" swivel casters.

Model No.	Part No.	Dimensions	Capacity	Ship Wt. (lb) est.
8766SS	0828766000	16"D 7-1/4"H	8.5 qt.	10


HAMPERS

Model No.	Part No.	Description	Ship Wt. (lb) est.
CHROME:			
2010	0962010000	Step-On Hamper, white poly lid	25*
2012	0962012000	Hamper Bag Storage Pouch	
2018	0962018000	Round 18" Hamper	11*
2025	0962025000	Round 25" Hamper	12*
2118	0962118000	Tilt-Top Hamper	21*
2230	0962230000	Folding Adjustable Hamper	18*

STAINLESS STEEL:

2301SS	0922301000	Triangular Space-Saver Hamper 18" D.	25*
7773SS	0927773000	Porter 26" diameter	20
7773SS-LP	0927773101	7773SS with plastic lid	35
7773SS-LFP	0927773201	7773SS foot-operated with plastic lid	55
7774SS	0927774000	Porter 18-5/8" diameter	11
7774SS-LF	0927774201	Porter 19-5/8" diameter	21
HB18	9839777201	18" diameter hamper bags	
HB44	9839775202	44" diameter hamper bags	
PLK-26	0927773548	Plastic Lid Kit-26"	

RECEPTACLES

Model No.	Part No.	Description	Ship Wt. (lb) est.
7766SS	0817766000	Lenox Kickbucket with 13 qt. pail	12*
7766SS-NP	0817766001	Lenox Kickbucket frame only	8*
KB-1	9708776600	13 qt. pail fits all models	4*
8766SS	0828766000	Round Sponge Receptacle	10*

* Items can be shipped via FedEx® Ground delivery service.



39 Robinson Road Lodi, NJ 07644 • 800-247-5070 • Fax: 201-909-0832
e-mail: info@blickman.com • www.blickman.com

All shipping weights are estimates only and may be changed without prior notice. Specifications may be improved and dimensions changed without prior notice. Unless otherwise noted, all products displayed or listed are deemed, 'MR-'Unsafe,' in accordance with ASTM Standard F2052-06.

FedEx service marks used by permission.

Copyright © 2006 Blickman, Inc.