

**BLICKMAN  
BUILT™**

# SURGICAL UTENSILS

Engineered for a Lifetime — Our Reputation is Stainless



DRESSING JARS  
EMESIS BASINS  
IODINE CUPS  
MEDICINE CUPS

OINTMENT JARS  
SOLUTION BASINS  
BUCKET PAIL  
CATHETER TRAYS  
GRADUATED MEASURE  
INSTRUMENT TRAY




**PRODUCT LISTING**
**SURGICAL UTENSILS**

Model No.	Part No.	Description	Capacity	Dimensions	Ship Wt. (lb)
DJ-1	9788000001	Dressing Jar	1/2 qt. (.5 liter)	4-1/8"C 2-5/8"D	1*
DJ-2	9788010002	Dressing Jar	1 qt. (1 liter)	4-1/8"C 5-1/4"D	1*
DJ-3	9788020003	Dressing Jar	2-1/4 qt. (1.9 liter)	4-7/8"C 6-7/8"D	2*
DJ-4	9788030004	Dressing Jar	3 qt. (2.8 liter)	6-1/8"C 7-1/4"D	2*
DJ-5	9788040005	Dressing Jar	4-1/2 qt. (4 liter)	7"C 7-3/8"D	2*
DJ-6	9788060006	Dressing Jar	6 qt. (5.7 liter)	7-1/2"C 9"D	2*
DJ-7	9788070007	Dressing Jar	8 qt. (7.6 liter)	8-1/4"C 10"D	3*
EB-1	9719090400	Emesis Basin	12 oz. (300cc)	6-5/8"W 3-1/2"D 1-1/2"H	1*
EB-2	9720595450	Emesis Basin	16 oz. (475cc)	8-1/8"W 3-5/8"D 1-3/4"H	1*
EB-3	9725511555	Emesis Basin	28 oz. (750cc)	10"W 4-1/2"D 2-1/8"H	1*
EB-4	9730514065	Emesis Basin	56 oz. (1500cc)	12"W 5-1/2"D 2-1/2"H	1*
IC-1	9775500002	Iodine Cup	6 oz. (150cc)	2"C 3-1/4"D	1*
IC-2	9710565001	Iodine Cup	14 oz. (355cc)	2-5/8"C 4-3/8"D	1*
MC-1	9784920001	Medicine Cup	2 oz. (50cc)	1-3/4"C 2"D	1*
OJ-1	9787720001	Ointment Jar	8 oz. (200 cc)	2-1/2"C 3-1/8"D	1*
SB-1	9737370120	Solution Basin	8-1/2 qt. (8 liter)	4-1/2"C 14-1/2"D	3*
1830	9701830000	Solution Basin	7 qt. (6.5 liter)	4-5/8"C 13-5/8"D	2*
KB-1	9708776600	Bucket Pail	13 qt. (12.3 liter)	10"C 11-1/2"D	4*
CT-1	9785030101	Catheter Tray w/ Cover		8-1/2"W 3"D 1-1/2"H	1*
CT-1T	9782831000	Catheter Tray (tray only)		8-1/2"W 3"D 1"H	1*
CT-1C	9782832000	Catheter Tray (cover only)		8-1/2"W 3"D 1/2"H	1*
GM-1	9795160001	Graduated Measure	16 oz. (475cc)		1*
GM-2	9795320002	Graduated Measure	32 oz. ( 950cc)		1*
GM-3	9795640003	Graduated Measure	64 oz. (1900cc)		1*
IT-1	9710015001	Instrument Tray		15"W 10"D	3*
IT-2	9713751912	Instrument Tray		9-1/8"W 13-3/4"D	4*
IT-3	9716002100	Instrument Tray		21"W 16"D	4*
IT-5	9736200250	Instrument Tray		25"W 20"D	6*
1505	9736150501	Instrument Tray		19-1/8"W 12-5/8"D	4*

\* Items can be shipped via FedEx® Ground delivery service.



39 Robinson Road Lodi, NJ 07644 • 800-247-5070 • Fax: 201-909-0832  
 e-mail: [info@blickman.com](mailto:info@blickman.com) • [www.blickman.com](http://www.blickman.com)