



Ambulatory Surgery Centers (ASC)

Quick, mobile, sterile environment

Includes:

- Cabinets and Casework
- Case Carts
- Chairs
- **Exam Stools**

IV Stands

- **Foot Stools**
- Hampers

- Mayo Stands
- Scrub Sinks
- **Solution Stands**
- **Surgical Utensils**
- Tables
- **Utility Carts**
- Warming Cabinets



Customizable, secure, and organized

Includes:

- Freestanding Cabinets
- Customizable Casework
- Sinks and Countertops
- Console Cabinets
- Accessories and Options




THE BLICKMAN **DIFFERENCE**

Freestanding Cabinets

- Our stainless steel cabinets are constructed with a seamless welded face preventing dirt collection while providing long lasting durability through unparalleled strength and corrosion resistance.
- Our freestanding cabinets rest on four legs with glides or can be mobilized on casters.
- Add a sloped top to prevent dirt collection to create a sterile environment.
- Kay and Kennedy cabinets feature two separate sections with unlocked or locked compartments.
- Paul and Walter cabinets provide maximum storage room for large equipment.

Key Product Features



Cabinets

Freestanding

FREESTANDING

Description: Our Kay and Kennedy freestanding cabinets feature separate upper and lower sections. The upper sections of these cabinets are locked, transparent compartments for easily viewing and securing contents. The Kennedy cabinet includes an independently locked and welded narcotic locker.



Kay Cabinet

Upper section: Consists of a locked 29 $^3/_4$ " x 31 $^3/_8$ " x 15" compartment with a glass door and three removable glass shelves

Lower section: Consists of a 29 $^3/_4$ " x 20 $^1/_2$ " x 15" compartment with a stainless steel door and one stainless steel shelf



Kennedy Cabinet

Upper Section: 8" deep and contains a welded and independently locked narcotic locker

Lower Section: 16" deep and contains three 7" deep drawers and one 15" deep compartment

Name	Kay	Kennedy
Dimensions (WxHxD)	30" x 58" x 16"	36" x 76 ½" x 16"
Model #	7953SS	7956SS
Item #	1517953000	1527956000
Shipping Weight	219 lbs.	291 lbs.

Options: 2 "Casters, stainless steel shelves for upper cabinets, sloping top, all stainless steel doors

Custom Designs Available Upon Request

Cabinets

Freestanding

FREESTANDING

Description: Our Paul and Walter cabinets are double door freestanding cabinets. They each come standard with two glass locked doors for securing contents. Our Walter cabinet is ideal for storing large instruments on its back pegboard, while our Paul cabinet includes shelving for storing lots of instruments and supplies. They are both 16" deep and 76" tall on gliders. They are 79 \frac{1}{4"} tall on casters. Choose from two standard width sizes: 35 \frac{5}{8"} and 47 \frac{5}{8"}.



Walter Cabinet

This cabinet is great for storing endoscopy equipment or other large instruments on its pegboard



Paul Cabinet

This cabinet provides maximum shelving storage for instruments and supplies. Glass shelves are adjustable every 1/2''

Name	Wo	Walter		ıul
Width	35 ⁵ /8″W	47 ⁵ /8″W	35 ⁵ /8″W	47 ⁵ /8"W
Model #	7959SS	7960SS	7970SS	797ISS
Item #	1537959000	1537960000	1537970000	1537971000
Shipping Weight	269 lbs.	310 lbs.	330 lbs.	400 lbs.

Options: Glides or casters, recessing, flat top, stainless steel shelves (Paul Cabinet), all stainless steel doors

Custom Designs Available Upon Request



FOOT OPTIONS



Locking Casters



Casters



Gliders

ACCESSORIES



Electronic Locking Systems, available for most of our cabinets



Additional Shelves, your choice of stainless steel or glass



Sloped top (see order details below)

Model#	Item#	Description	Shipping Weight
ST-3016	4503516300	Sloped Top for Kay (30" x 6" x 16")	35 lbs.
ST-3608	4503536080	Sloped Top for Kennedy (36" x 6" x 8")	37 lbs.
ST-3516	4503516350	Sloped Top for Walter or Paul (slim) (35 ⁵ /8" x 6" x 16")	43 lbs.
ST-4716	4503516470	Sloped Top for Walter or Paul (wide) (47 ⁵ /8" x 6" x 16")	54 lbs.



	Construction Specifications
Standard Dimensions	All Freestanding Cabinets: 16" deep
Materials	304 stainless steel ¹ / ₄ " clear tempered sheet glass
Gauge	Body, Sloped Top, Drawers, and Shelf Inserts: 18GA; Legs: 12GA Glass Doors (Front and Back): 18GA Solid Doors (Front): 18GA; Solid Doors (Back): 20GA
Welding	All welding by heliarc method and exposed welds are ground smooth
Finish	Exposed surfaces are free from pit marks, weld seams, scale and are polished with #4 polish to a satin finish
Casters	2" or 3" ball bearing, swivel casters with brakes
Notes	 Kay: A counter weight is fastened underneath at bottom rear of cabinet for doors in open position. Upper door fitted with lock handle and three way latching device engaging at top.
	• Kennedy: Right hand door equipped with a locking knob and three way latching. Left hand door with dummy knob and roller latch. An independently locked narcotic cabinet, 13 1/4" x 5" x 7 3/4" with double walled door is welded to the bottom, back, and side of the compartment.
	Paul & Walter: A counter weight is fastened underneath at bottom rear of cabinet for doors in open position. Right hand door is fitted with a locking handle and three point latching, and the left door has a dummy knob and spring loaded nylon roller catches.


Shipping Information

All shipping weights are estimates and may be changed without prior notice. Specifications may be improved and dimensions changed without prior notice. Please contact customer service for estimated shipping costs at info@blickman.com.

THE BLICKMAN **DIFFERENCE**

Casework

- Need more flexibility in today's space squeezed environment? We provide you with casework that can be assembled together and modified to fit your custom dimensions.
- Whether it is renovation or new construction, we will be happy to work with your facility's personnel, architects, contractors, and/or equipment planners to meet your custom requirements.
- Utilizing 3D Design Software REVIT®, we have developed a library of our standard products which is available at no charge to your company. Our 3D library will save you considerable time and money when working on projects that require stainless steel cabinets and equipment.
- Each product line of cabinets (console cabinets, high cabinets, base cabinets, and wall cabinets) have their individual properties identified for easy customization: (cubic feet, weight, "X", "Y" and "Z" dimensions, voltage, etc.).






HIGH CABINETS

Description: Our stainless steel high cabinets are designed to stand alone or fit within a framework of cabinets. They are constructed as floor and wall supported individual cases with finished sides. They each contain five adjustable stainless steel shelves, with the option of glass shelves and additional shelves. They are all 84" tall and 18" deep cabinets. Width varies and other alternative dimensions are available.

Customize the widths and depths of these models to your room.



Name	Single Door	Double Door	Double Sliding Doors
		Glass Doors	
Dimensions (WxHxD)	18 ³ / ₄ " x 84" x 18"	35" x 84" x 18"	35" x 84" x 18"
Model #	K18HAS	L35HAS	M35HAS
Item #	2030518000	203053500	2030635000
Shipping Weight	215 lbs.	310 lbs.	325 lbs.
		Solid Stainless Steel Doors	3
Dimensions (WxHxD)	18 ³ / ₄ " x 84" x 18"	35" x 84" x 18"	35" x 84" x 18"
Model #	H18HAS	I35HAS	J35HAS
Item #	2030418000	2030435000	2030535180
Shipping Weight	230 lbs.	300 lbs.	347 lbs.

Options: Sloped top, glass shelves, additional shelves

Custom Designs Available Upon Request




Customizable Casework

WALL CABINETS

Description: Our stainless steel wall cabinets are designed to stand alone with finished sides or to fit within a framework of cabinets. They are constructed for flush face wall surfaces with a recessed back panel for uneven wall surfaces. Each cabinet comes with two 1/2" adjustable stainless steel shelves. Choose from swinging or sliding doors made of either solid stainless steel or glass. Our wall cabinets are all 30" tall and 13" deep. Width varies and other alternative dimensions are available.

Customize the widths and heights of these models to your room.

Common Dimensions:



Custom Designs Available Upon Request

Cabinets

Customizable Casework

OPEN WALL SHELVES

Description: Our stainless steel open wall cabinet is designed to stand alone with finished sides or fit within a framework of cabinets and comes with two adjustable stainless steel shelves. Alternatively, we offer a wall shelf that can be placed anywhere along a wall.



Common Dimensions:

Wall Cabinet Widths: All Single nd Double Width Sizes Heights: 36" and 48" Wall Shelf Widths: 48", 60", and 72", Depths: 16"



Name Width	Open Wall Cabinet		Wall Shelf	
	35″ W	47" W	36" W	84" W
Dimensions (WxHxD)	35" x 30" x 13"	47" x 30" x 13"	36" x 12"	84" x 12"
Model #	G35LS	G47LS	WS3612	WS8412
Item #	2020435000	202044700	1813612000	1818412000
Shipping Weight	75 lbs.	142 lbs.	26 lbs.	40 lbs.

OPEN BASE CABINETS

Description: Our stainless steel open base cabinets offer open shelving to fit within a framework of cabinets or stand alone with finished sides. They are constructed for flush face wall surfaces and come with two adjustable stainless steel shelves. Our open base cabinets are all 22" deep and 35 3/4" high and come in a variety of width sizes.

Other Commonly Ordered Single Width Sizes: 21 1/8" W
Other Commonly Ordered Double Width Sizes: 29"W and 41"W





Name	Single Widtl


Width	18 ¹/ ₈ " W	24 ¹/8″ W	35″ W	47" W
Model #	AJ18HS	AJ24HS	AJ35HS	AJ47HS
Item #	2012518000	2012524000	2012535000	2012547000
Shipping Weight	102 lbs.	122 lbs.	152 lbs.	170 lbs.

Custom Designs Available Upon Request



SINGLE BASE CABINETS

Description: Our stainless steel single width base cabinets come in five different configurations and are designed to stand alone with finished sides or to fit within a framework of cabinets. They are all 35 3/4" tall and 22" deep, and come in two width sizes, 18 1/8" and 24 1/8". **Other Commonly Ordered Single Width Sizes:** 21-1/8"W



Note: Cabinets without vents include one stainless steel shelf inside



183 lbs.

Name	-AS-	-AE-	
	3 Drawer	4 Drawer	
	18 1/8	₃ " Wide	
Model #	AS18HS	AE18HS	
Item #	2012718000	2012318000	
Shipping Weight	167 lbs.	174 lbs.	
	24 1/:	₈ ″ Wide	
Model #	AS24HS	AE24HS	
Item #	2012724000	2012324000	

Note: Cabinets without vents include one stainless steel shelf inside

Shipping Weight

192 lbs.

Cabinets

Customizable Casework

DOUBLE BASE COMPARTMENT CABINETS

Description: Our stainless steel double base compartment cabinets come in seven different configurations and are designed to stand alone with finished sides or to fit within a framework of cabinets. They are all 35 3/4" tall and 22" deep, and come in two width sizes, 35" and 47". **Other Commonly Ordered Sizes:** 29"W and 41"W









Name	-DA- Double Door	-DB- Double Vent (Sink)	-DC- Large Drawer	-DD- Double Drawer
		35″ V	Vide	
Model #	DA35HS	DB35HS	DC35HS	DD35HS
Item #	2013235200	2013235000	2013235001	2013335000
Shipping Weight	161 lbs.	145 lbs.	200 lbs.	194 lbs.
		47″ V	Vide	
Model #	DA47HS	DB47HS	DC47HS	DD47HS
Item #	2013247200	2013247000	2013247001	2013347000
Shippina Weight	205 lbs.	178 lbs.	235 lbs.	243 lbs.

Note: Cabinets without vents include one stainless steel shelf inside







-DG-	-DJ-	-DK-
1 Door, 4 Drawer	1 Door, 8 Drawer	1 Door, 5 Drawer
	35" Wide	
DG35HS	DJ35HS	DK35HS
2013435000	2013535001	2013535000
180 lbs.	225 lbs.	257 lbs.
	47" Wide	
DG47HS	DJ47HS	DK47HS
2013447000	2013547001	2013547000
282 lbs.	280 lbs.	283 lbs.
	DG35HS 2013435000 180 lbs. DG47HS 2013447000	35" Wide DG35HS DJ35HS 2013435000 2013535001 180 lbs. 225 lbs. 47" Wide DG47HS DJ47HS 2013447000 2013547001

Note: Cabinets without vents include one stainless steel shelf inside



Customizable Casework

DOUBLE BASE DRAWER CABINETS

Description: Our stainless steel double base drawer cabinets come in five different configurations and are designed to stand alone with finished sides or to fit within a framework of cabinets. They are all 35 3/4" tall and 22" deep, and come in two width sizes, 35" and 47".

Other Commonly Ordered Sizes: 29"W and 41"W





Name

Shipping Weight

-DP-4 Drawer

276 lbs.

5 Drawer

315 lbs.

	35	" Wide	
Model #	DP35HS	DS35HS	
Item #	2013735000	2013735001	
Shipping Weight	238 lbs.	275 lbs.	
	47	" Wide	
Model #	DP47HS	DS47HS	
Item #	2013647000	2013747001	







Name	-DT-	-DR-	-FC-
	8 Drawer	4 Small, 2 Large	6 Small, 2 Large
		35" Wide	
Model #	DT35HS	DR35HS	FC35HS
Item #	2013835000	2013635000	2013235002
Shipping Weight	287 lbs.	265 lbs.	275 lbs.
		47" Wide	
Model #	DT47HS	DR47HS	FC47HS
Item #	2013847000	2013747000	2013247002
Shipping Weight	324 lbs.	310 lbs.	323 lbs.

Custom Designs Available Upon Request




SINKS AND COUNTERTOPS

Please email us at info@blickman.com



Description: Our sinks and countertops must be ordered separately and are installed in the field. Most of our countertops are fabricated from one piece of fully welded, 16 gauge, 304 stainless steel. They extend 1" beyond the front, right, and left sides of base cabinet dimensions. Unless otherwise specified, they come with a Marine Edge to help protect against liquid slippage. Optionally, 4" high backsplashes and endsplashes are available to protect wall surfaces.

Base Cabinets Compatible with Additional Sink Bowl




We suggest using Models AC and DB for moisture ventilation.

All piping and sink bowls can fit under any of these models with the door compartments and no drawers on top.




DRAWER HOUSING

Please email us at info@blickman.com



Description: You can use drawer housing to add a drawer or two under your countertop or tabletop surface. This helps increase storage space to prevent lack of access to equipment and supplies when stored separately. All drawer housing is 22" deep and reaches 5" below the underside of the countertop or tabletop. Width of drawer housing varies depending on size and number of drawers inside. The drawers are all 18 1/8" deep and 2 1/2" tall and come in two width sizes: 14 1/8" and 20 1/8".





Name	Single	Drawer	Double	Drawer
Drawer & Housing Widths	14 ¹ / ₈ ″W Drawer 24″W Housing	20 ¹ / ₈ "W Drawer 29"W Housing	14 ¹ / ₈ "W Drawer 35"W Housing	20 ¹ / ₈ "W Drawer 47"W Housing
Model #	DH124S	DH129S	DH235S	DH247S
Item #	2043424000	2043429000	2043435000	2043447000
Shipping Weight	35 lbs.	40 lbs.	50 lbs.	60 lbs.

Custom Designs Available Upon Request

OPTIONS



Sloped Top is available for all High Cabinets and Wall Cabinets. Listed below are the sloped tops that correspond to the cabinets with the standard dimensions shown in this previously in this section.

Model#	Item#	Description	Shipping Weight
ST-1818	5003518180	Sloped Top for 18 ³ / ₄ "W High Cabinets (18 ³ / ₄ " x 11 ¹ / ₄ " x 18")	24 lbs.
ST-3518	5003535180	Sloped Top for 35"W High Cabinets (35" x 11 1/4" x 18")	41 lbs.
ST-1813	5003518130	Sloped Top for 18 ³ / ₄ " W Wall Cabinets (18 ³ / ₄ " x 8 ¹ / ₄ " x 13")	21 lbs.
ST-3513	5003535130	Sloped Top for 35" W Wall Cabinets (35" x 8 1/4" x 13")	36 lbs.

ACCESSORIES



Electronic Locking Systems available for most of our cabinets

Additional Shelves, your choice of stainless steel or glass

Custom Designs Available Upon Request



	Construction Specifications
Standard Dimensions	High Cabinets: 84" tall x 18" deep, widths vary Wall Cabinets: 30" tall x 13" deep, widths vary Base Cabinets: 35 3/4" tall x 22" deep, widths vary Countertops: 1" overhang in front and both sides of base cabinets
Materials	304 Stainless Steel 1/4" clear tempered sheet glass
Gauge	Body, Sloped Top, Drawers, and Shelf Inserts: 18GA; Base: 16GA Glass Doors (Front and Back): 18GA Solid Doors (Front): 18GA; Solid Doors (Back): 20GA
Welding	All welding by heliarc method and exposed welds are ground smooth
Finish	Exposed surfaces are free from pit marks, weld seams, scale and are polished with #4 polish to a satin finish
Casters	N/A
Notes	 Cabinet Parts: Hinge Doors: Hung on semi concealed stainless steel hinges and fitted with die cast aluminum bar pulls and spring loaded nylon roller catches. Sliding Doors: Suspended from overhead ball bearing rollers, which are noiseless and self-lubricating, and ride in stainless steel guide strips fitted with die formed, recessed pulls securely welded in place. Shelves: Constructed with a channel shaped, flanged edge and centers the full height of the cabinet. Drawers: Drawer bodies are a one piece construction with covered corners. Drawer face equipped with die cast aluminum bar pole handles. Drawers operated on ball bearing roller slides with full extension. High Cabinets & Wall Cabinets: Recessed pan back, 3/4" deep. Base Cabinets: Each applied with a channel shaped base, incorporating toe space approximating 3" deep and 4" high, for easy cabinet access. Can be omitted when recessed or another base is used. Sinks and Countertops: Fabricated in one piece, where possible, with welded corners, channel shaped exposed edges, and reinforced with three full length hat shaped galvanneal channels.

All shipping weights are estimates and may be changed without prior notice. Specifications may be improved and dimensions changed without prior notice. Please contact customer service for estimated shipping costs at info@blickman.com.


Shipping Information

THE BLICKMAN **DIFFERENCE**

Console Cabinets

- Our console cabinets are designed with you in mind. Without taking up valuable floor space, choose these cabinets for storing a full scheduled day, protecting your sterile supplies behind closed doors, and allowing instant access to supplies and equipment to minimize procedural delays.
- We offer individual units that can be combined and assembled into your own custom console arrangements. You can easily design your own by selecting cabinets with solid or glass doors available with or without locks. Locks can be keyed alike or independently. Stainless steel perimeter recessing trim may be purchased with all built-in installations.
- We offer units with shelving, good for storing equipment and supplies and other units with pegboard backs, great for storing large pieces of equipment.
- We also offer a desk cabinet with drawers, a sloped writing surface, and a dry erase board.

Key Product Features



Wide Shelving Cabinet provides maximum shelving storage for instruments and supplies


Doors can be a transparent glass or solid stainless steel




Console Cabinets

SINGLE WIDTH

Description: Our single width console cabinets are 24 1/8" x 60" x 18" and are designed as individual units that can be combined and assembled into your own custom console arrangements. All doors and shelves come in glass or solid stainless steel. Drawers are all stainless steel.





Name	Miscellaneous Supply Cabinet		Desk Console Cabinet
	Glass Door	Solid Door	Drawers
Model #	MS24	MS24S	DC24
Item #	1716724000	1716724258	1713224000
Shipping Weight	240 lbs.	240 lbs.	245 lbs.

Custom Designs Available Upon Request



Console Cabinets

DOUBLE WIDTH

Description: Our double width console cabinets are all designed as individual units that can be combined and assembled into your own custom console arrangements. All doors and shelves come in glass or solid stainless steel. They are both 60" high and 18" deep and come in two width sizes: 35" and 47".





name	Equipment & Suppl	ies Console Cabinet	lable Accessorie	s Console Cabinet
	Glass Door	Solid Door	Glass Door	Solid Door
		35″ \	N	
Model #	ES35	ES35S	TA35G	TA35
Item #	1713735000	1713735258	1718235010	1718235000
Shipping Weight	250 lbs.	250 lbs.	250 lbs.	250 lbs.
		47" \	N	
Model #	ES47	ES47S	TA47G	TA47
Item #	1713747000	1713747258	1718247003	1718247000
Shipping Weight	275 lbs.	275 lbs.	355 lbs.	355 lbs.

Custom Designs Available Upon Request

Cabinets

Console Cabinets

RECESSED WARMING CABINETS

Description: Our recessed warming cabinets are designed for built-in installation. Single chamber glass door, solid door, and passthrough models available. Three perforated stainless steel shelves come standard in these models.

Glass Door



Solid Door



Passthrough



Note: For the passthrough unit, the front door is glass, the back door is solid.

Dimensions

Width: 30" x Height: 60" x Depth: 20 5/8"

Electrical:

Amps: 7.5 | Total Wattage: 750

Name	Glass Door	Solid Door	Passthrough
Model #	SW40TG	SW40TS	SW40TPTGS
Item #	14BSW40243	14BSW40200	14BSW40278
Shipping Weight	320 lbs.	320 lbs.	245 lbs.

Options: Dual chambers, alternative passthrough configurations.



ACCESSORIES



Recess trim available for built-in console cabinets

Model#	Description	Shipping Weight
CG4448	48" High corner guards, for protection and corrosion resistance (4"x4"x48")	20 lbs.
CG4496	96" High Corner guards, for protection and corrosion resistance (4"x4"x96")	20 lbs.
Rt120	Recessing trim, 22GA 304 stainless steel (10' strips, $1^{1}/_{2}$ " x $1^{1}/_{2}$ wide)	20 lbs.



Electronic Locking Systems available for most of our cabinets



Passthrough Cabinets available for built-in applications

Custom Designs Available Upon Request



	Construction Specifications
Standard Dimensions	All console cabinets are 16" high x 18" deep, widths vary Exception: warming cabinets are 20" deep
Materials	304 Stainless Steel 1/4" clear tempered sheet glass
Gauge	Body, Sloped Top, Drawers, and Shelf Inserts: 18GA Glass Doors: Front and Back: 18GA Solid Doors (Front): 18GA; Solid Doors (Back): 20GA
Welding	All welding by heliarc method and exposed welds are ground smooth
Finish	Exposed surfaces are free from pit marks, weld seams, scale and are polished with #4 polish to a satin finish
Casters	N/A
Notes	 All Console Cabinets: Solid doors are double panned, channel reinforced. Glass doors come in double shell design, with one piece outer frames dieformed to receive rubber gasketed glass panels held in place by inner retaining strips. Double Door Console Cabinets: Right door fitted with two point up and down knob latching, left hand door with integral stops, dummy knob, and spring loaded nylon roller catch. Doors operate on three semi-concealed stainless steel hinges. Pegboard Cabinet interiors (back and sides) are lined with removable 18GA stainless steel pegboard. Desk Console Cabinet: Four drawers with double wall fronts operate on side mounted, two section ball bearing roller slides with built-in extension stops. Heads provided with diecast aluminum bar pulls. Single Door Shelving Console Cabinet: Can be either right or left hinged. Recessed Warming Cabinets: Temperature control from 90°F - 160°F. 160×180 pixel graphical display with dual colors (red/green). Built-in temperature data recorder with real time clock (RTC).

Shipping Information

All shipping weights are estimates and may be changed without prior notice. Specifications may be improved and dimensions changed without prior notice. Please contact customer service for estimated shipping costs at info@blickman.com.

Blickman REVIT® Library

We have developed a library of our standard products using 3D Design Software REVIT®, which is available to your company at no charge.

Included in our Revit Library you can find:
Free standing Cabinets
Casework
Case Carts
Tables
Warming Cabinets

Visit: www.blickman.com/revit-files



Create your entire room from our cabinet files to discover a design that works best for your needs. We are happy to work with you to create custom cabinetry, casework, tables, etc. to fit in your space.

Custom Designs Available Upon Request



Case Carts

Compact, stable, secure transport

Includes:

- Closed Case Carts
- Open Case Carts
- Shelving Options
- Accessories



THE BLICKMAN **DIFFERENCE**

- Operating room surgical staff can transport fully loaded case carts throughout the facility with the doors closed for increased protection and cart cleaning is a simple and easy task.
- We offer options that meet a variety of different surgery needs including neurosurgery, general surgery, outpatient surgery, labor and delivery, and more. We also offer space saving solutions that take up less floor space while being able to handle more extensive procedures.
- Our broad offering of case carts remains the most durable and highest functioning on the market. We are the only manufacturer in our class to offer a lifetime limited warranty.
 Our advanced engineering capabilities enable us to customize case cart designs to satisfy virtually any requirement your medical facility may have.
- The doors of the closed case carts swing a full 270° and a positive spring acting clip at the cart ends capture doors to stay in the open position during cart washing and while in use in the operating room.

Key Product Features





Maxi to Mini

Description: Our while maxi, multipurpose, and mini case carts differ in width, they are all 40 1/2" tall and 29" deep.. They each come standard with one wire shelf. The maxi case cart comes standard with a heavy duty roll-out shelf, and the multipurpose and mini case carts come with a lightweight pull-out shelf.

Maxi Cart

52"W

Ideal for procedures that require a larger volume of supplies such as neurosurgery and orthopedic surgery

Internal Dimensions: 42 7/8"W x 29 3/8"H x 25"D

Multipurpose Cart

42" W

Ideal size for most routine surgical procedures, whether or not your facility uses a container system

Internal Dimensions: 32 7/8"W x 29 5/8"H x 25"D

Mini Cart

29 5/8" W

Designed for labor and delivery and outpatient surgery, where a smaller cart is preferred

Internal Dimensions: 20 1/2"W x 29 5/8"H x 25"D







1	
	6

Push Handle



Item #

Floor Space

Shipping Weight

Name	Maxi	Multipurpose	Mini

CCC2-19	CCC3-19			
CCC2-19	CCC3-19			
2293332000	2293333000			
8 ¹ / ₂ ft ²	6 ft²			
225 lbs.	168 lbs.			
Glass Doors				
0000 100	CCC3-19G			
_	CCC2-19G			

22933G2000

8.46 ft²

219 lbs.

Custom Designs Available Upon Request

22933G1000

10.47 ft²

260 lbs.

Please email us at info@blickman.com

22933G3000

5.97 ft²

165 lbs.

Case Carts

Space Saving

Description: Our space saver case carts are designed to give a greater capacity while using the same floor space as the mini and multipurpose case carts. Their height and additional shelving provide additional storage capacity. Instead of 40 1/2" tall, they are 55 1/2" tall, and instead of one shelf, they come standard with two lightweight pull-out wire shelves.

Space Saving Cart

29 ⁵/₈" W

This cart can handle most routine surgical procedures.

Recessed Handle





Name Space Saving

Ultra Space Saving Cart 42" W

This cart can handle more extensive procedures.

Internal Dimensions: 32 7/8"W x 44 5/8"H x25"D



Ultra Space Saving





	Solid	Doors
Model #	CCC4-19	CCC5-19
Item #	2293334000	2293335000
Floor Space	6 ft²	8 ¹ / ₂ ft²
hipping Weight	208 lbs.	276 lbs.
	Glass	Doors
Model #	CCC4-19G CCC5-19G	
Item #	22933G4000	22933G5000
Floor Space	6 ft²	$8^{1/2} ft^{2}$
nipping Weight	198 lbs.	262 lbs.

Custom Designs Available Upon Request

Case Carts

Over The Road

Description: As healthcare processes change, Blickman works with customers to design the products they need. With growing demand for off-site sterile processing, Blickman reengineered its case cart lineup by adding a new series of carts dedicated to meet its customers' new and more demanding transport requirements.



*Ultra Space Saver - OTR5 is shown here with optional vertical side handles.



Name	Item#	Model#	Dimensions (WxHxD)	Floor Space	Shipping Weight
Maxi	2327331030	OTRI	53 x 42 x 29	10 1/2 ft²	266 lbs.
Multipurpose	2327332030	OTR2	43 x 42 x 29	8 1/2 ft ²	225 lbs.
Mini	2327333030	OTR3	30 5/8 x 42 x 29	6 ft²	168 lbs.
Space Saver	2327334030	OTR4	30 5/8 x 56 3/4 x 29	6 ft²	208 lbs.
Ultra Space Saver	2327335030	OTR5	43 x 56 3/4 x 29	8 1/2 ft²	276 lbs.

Custom Designs Available Upon Request



Open

Description: Our open case cart designs are easy to clean, dry, load, and maneuver. The full encircling bumper and solid bottom prevents cross contamination. This design is ideal for carts that go directly to a sterile core, bypassing non-sterile transport.

Overall Dimensions:

Width: 42" x Height: 42" x Depth: 29"





Name	Solid Shelves	Wire Shelves	
Model #	OCC3	OCC4	
Item #	2410396000	2410496000	
Shipping Weight	219 lbs.	141 lbs.	

Custom Designs Available Upon Request

Case Carts

Interior Shelving

STANDARD

Wire Roll Out Shelf



Wire Pull Out Shelf



CCC1 CCC2, 2E, 3, 4, and 5

OPTIONAL





Roll Out Shelving

Pull Out Shelving

Case Carts	Wire	Perforated	Solid	Wire	Perforated	Solid
CCCI	611C1RWKIT	611C1RPKIT	613C1RSKIT	N/A	61517P1KIT	61517S2KIT
CCC2/CCC5	610C2RWKIT	610C2RPKIT	610C2RSKIT	6152272KIT	6152275KIT	6152274KIT
CCC3/CCC4	610C1RWKIT	610C1RPKIT	610C1RSKIT	6152273KIT	6152271KIT	6152073KIT



Exterior Shelving

EXTENSION SHELVES



Available on Model CCC2E Only

- Expand your work surface area -

DROP SHELVES



Compantible with CCC1, CCC2, & CCC3

- Expand your work surface area -

Custom Designs Available Upon Request

Case Carts

Accessories

ACCESSORIES

Clean/Dirty Indicators







US Design Patent #: US D883,378 S US Utility Patent#: US 10,643,505 B1

US Design Patent #: US D883,379 S US Utility Patent#: US 10,643,505 B1

Part#	Description	Shipping Weight	
5278628KIT	Angle Model	2 lbs.	
5200045KIT	Dial Model	2 lbs.	
5278618KIT	Pod Model	2 lbs.	





Card Holder

We offer a 8.5" x 11" card holder. Custom sizes and placement options available. Stainless steel frame withstands prolonged use and cutouts allow for drainage.

Tote Boxes

Sturdy, long lasting tote boxes are ideal for organizing and storing small supplies that may slip through the spaces of wire shelves.

Part#	Description	Shipping Weight	
5206701185	Card Holder	N/A - Factory Installed	
9854194720	Tote Boxes	5 lbs.	

Construction Specifications

Universal Dimensions	Closed Case Carts: 29" deep, varying widths. Standard Height: 40 1/2"; Space saving height: 55 1/2" Open Case Carts: 42" x 42" x 29"
Materials	304 stainless steel, non-magnetic
Gauge	Closed Case Carts: Body & Doors: 18GA; Top & Chassis: 16GA Open Case Carts: Top & Bottom: 16GA; Shelf: 18GA Bumper Channel: 14GA
Welding	All welding by heliarc method and exposed welds are ground smooth
Finish	Exposed surfaces are free from pit marks, weld seams, scale and are polished with #4 polish to a satin finish
Casters	Four 6" stainless steel swivel casters with polyurethane wheels, weight rated at 300lbs (one with directional lock)
Notes	Found on next page

Shipping Information

All shipping weights are estimates and may be changed without prior notice. Specifications may be improved and dimensions changed without prior notice. Please contact customer service for estimated shipping costs at info@blickman.com.

Custom Designs Available Upon Request

Notes Continued

CLOSED CASE CARTS:

- **Body:** One piece top and one piece body. Sound deadening angle welded under the top. Sound deadening channel on the center back. A pan shaped, step down bottom with a slope that facilitates draining. Horizontal intersections filled with heat resistant silicone. The body shell is bolted to the chassis. A tubular push handle is provided at one end of the cart.
- **Doors:** Double walled, channel enforced, (double or triple) hinged doors with paddle latch pulls on leading edges. Doors less than 36" high swing on two hinges, while doors greater than 36" high swing on three hinges. All doors permit full 270° swing. In fully opened position, positive spring acting clips at the cart ends capture doors for cart washing or operating room use.
- **Shelves:** Shelves slide on adjustable stainless steel wires, suspended from keyhole punched standards welded to the cart interior.

OPEN CASE CARTS:

- **Frame:** U-shaped tubular framing formed around 3" radii or tubular corner posts stabilized and capped by a pan shaped top welded to posts at contact points. Posts are fitted with casters and either rolling bumpers or encircling rubber bumpers.
- Tops, bottoms, fixed shelves: Fabricated with turned down edges, either hemmed or flanged, with corners notched to sit over or snuggly against the inside contour of the upright members. Notched corners are fully welded to the uprights the full height of the facing edges as well as around the notch.
- **Adjustable shelves:** Solid rods welded between the tubular uprights at each end, stabilizing the framing members while providing varying levels of shelf adjustment.

Custom Designs Available Upon Request



Chairs

Providing ultimate comfort and function

Includes:

- Patient Room Chairs
- Waiting Room Chairs



THE BLICKMAN DIFFERENCE

Chairs

- Our patient chairs come in our same washable vinyl (75 colors available see upholstery catalog) and are designed to be used in patient rooms and offices that require easy-to-clean, sterile environments.
- Our waiting room chairs come in our new 100% Polyester fabric (19 colors available see upholstery catalog) and are designed to be used in more comfortable and inviting environment such as waiting areas.
- All of our chairs include a backrest and four black washable plastic legs. We offer a standard size with and without arm rests.

Key Product Features





Patient & Waiting Room

Patient Room Chairs: Our patient room chairs come in our vinyl fabric for comfortable seating designed for use in patient exam rooms, offices, hallways, and other areas of your facilities that require sterilization.

Waiting Room Chairs: Our waiting room chairs come in our polyester fabric for more comfortable and inviting environments such as waiting areas.

Sizes: Our standard chair size is 24" x 33" x 22 5/8", comes with or without arm rests and has a weight capacity of 300 lbs.





Name	Chair	Chair with Arms	
Fabric	Vinyl (Patient Room Chairs)		
Model #	1130	1130WA	
Item #	1051130000	1051130100	
Shipping Weight	18 lbs.	21 lbs.	
Fabric	Polyester Fabric (Waiting Room Chairs)		
Model #	1120	1120WA	
Item #	1051120000	1051120100	
Shipping Weight	18 lbs.	21 lbs.	

Options: 75 vinyl colors and 19 polyester fabric colors- see upholstry catalog



Construction Specifications

Standard Dimensions	Standard Seat: Width: 24" x Height: 33" x Depth: 22 5/8"
Materials	Arms and Legs: Black Plastic
	Seats: 75 vinyl colors or 19 polyester fabric colors - see upholstery catalog
Gauge	N/A
Welding	N/A
Finish	N/A
Casters	N/A

Shipping Information

All shipping weights are estimates and may be changed without prior notice. Specifications may be improved and dimensions changed without prior notice. Please contact customer service for estimated shipping costs at info@blickman.com.

Toll Free: (800) 247-5070 Local: (973) 330-0557



www.blickman.com info@blickman.com

BLICKMAN



Providing ultimate comfort and function

Includes:

- Spindle Adjusted Stools
- Pneumatic Stools
- Parts & Accessories



THE BLICKMAN DIFFERENCE

EXAM STOOLS

- Our durable, high quality chrome and stainless steel exam stools come with 3" thick, dense foam 15" diameter cushion for ultimate comfort.
- Our upholstery is washable vinyl and is available in an array of 75 fashion colors that meet specific design requirements see upholstery catalog for details.
- For applications that require the longevity that stainless steel provides, pick from our revolving stools, or for a more modern look, choose from four or five leg models with spindle adjusted or pneumatic lift mechanisms.
- For back support, choose one of our models with a back rest.
- Our three styles of locking casters allow for easy mobility and safety.

Key Product Features





Spindle

4 LEG SPINDLE EXAM STOOLS

Description: We offer two four leg spindle adjusted exam stools with the option of 2" casters or rubber crutch tips. Choose our Passaic stool and sit comfortably on a 15" diameter, 3" thick polyfoam padded seat upholstered in vinyl or choose our Clifton stool for an easy-to-clean stainless steel seat with a non-slip vinyl insert. Add a removable pad slipcover to your Clifton stool for a more comfortable option.



		Standard	With Pad	Standard	
_			Casters		
	Model #	7745SS-C	7745SS-PC	7714-PSSC	
	Item #	1027745001	1027745003	1027714125	
	Shipping Weight	17 lbs.	17 lbs.	17 lbs.	
			Crutch Tips		
	Model #	7745SS	7745SS-P	7714-PSS	
	Item #	1027745000	1027745002	1027714025	
	Shipping Weight	16 lbs.	16 lbs.	16 lbs.	

CLIFTON

Options: Available in standard black. Consult factory for requests in other colors for large order quantities

Name

PASSAIC



Spindle

5 LEG SPINDLE EXAM STOOL

Description: Our five leg spindle exam stool rests on a 24" diameter base. Choose between a black composite base or an aluminum base with Wall Saver™ tips on the ends to protect against scuff marks and offer a textured surface for foot grip. Sit comfortably on its 15" diameter, vinyl upholstered seat with 3" thick polyfoam padding. This model is available with or without a back.



	Name	5 Leg Spindle Stool	5 Leg Spindle Stool with backrest
		Aluminun	n Base
5	Model #	1113	1113WB
	Titem #	1021113000	1021113200
9	Shipping Weight	22 lbs.	26 lbs.
		Black Compo	osite Base
X	Model #	1113BLK	1113WB-BLK
5	Item#	1021113100	1021113300
-	Shipping Weight	22 lbs.	26 lbs.

Options: Foot ring, 75 different vinyl colors - see upholstery catalog

Exam Stools

Pneumatic

FOOT OPERATED EXAM STOOL

Description: Our five-leg foot operated exam stool have a push ring under a 24" diameter aluminum base with Wall Saver™ tips on the ends to protect against scuff marks and offer a textured surface for foot grip. Sit comfortably on its 15" diameter, vinyl upholstered seat with 3" thick polyfoam padding. This model is available with or without a back.



Name	Foot Operated Stool	Foot Operated Stool with Backrest
	Aluminum Base	
Model #	1205	1206
Item#	1041205000	1041206000
Shipping Weight	22 lbs.	28 lbs.

Options: Foot ring, 75 different vinyl colors - see upholstery catalog



Pneumatic

HAND OPERATED EXAM STOOLS

Description: Our five-leg hand operated exam stool rests on a 24" diameter base. Choose between a black composite base or an aluminum base with Wall Saver™ tips on the ends to protect against scuff marks and offer a textured surface for foot grip. Sit comfortably on its 15" diameter, vinyl upholstered seat with 3" thick polyfoam padding. This model is available with or without a back.

Standard Stools



	Name	Hand Operated Stool	Hand Operated Stool with Backrest
		Alumin	um Base
5	Model #	1210	1210WB
	Item #	1041210000	1041210200
•	Shipping Weight	22 lbs.	26 lbs.
		Black Com	posite Base
	Model #	1210BLK	1210WB-BLK
5	Item #	1041210100	1041210300
•	Shipping Weight	22 lbs.	26 lbs.

Options: Foot ring, 75 different vinyl colors - see upholstery catalog



Exam Stools

Pneumatic

HAND OPERATED EXAM STOOLS

Economy Stools

Description: Our economy hand operated exam stool rests on a 24" diameter black composite base with black plastic twin wheel casters. Sit in comfort on a 15" diameter, vinyl upholstered seat with 3" thick polyfoam padding. this model is available with and without a back.

Adjusts 7" from 18 $^{3}/_{8}$ " to 25 $^{3}/_{8}$ " off the ground







Name E	conomy	Stoo
--------	--------	------

	Without Back	With Back
Model #	1200	1201
Item #	1041200100	1041201300
Shipping Weight	22 lbs.	26 lbs.

Options: Foot ring, 73 different vinyl colors - see upholstry catalog

Deluxe Stools

Description: Our deluxe hand operated exam stool comes standard with a foot ring for leg-joint comfort. These stools rest on a 24" diameter black composite or aluminum base, which comes with Wall Saver™ tips to protect against scuff marks and offer a textured surface for foot grip. Both bases travel easily on 2" hard rubber casters. Sit in comfort on a 15" diameter, vinyl upholstered seat with 3" thick polyfoam padding with a backrest.

Adjusts 7" from 21 $\frac{3}{8}$ " to 28 $\frac{3}{8}$ " off the ground







Name	Deluxe Stool		
	Aluminum Base	Black Composite Base	
Model #	1212	1212BLK	
Item #	1041212000	1041212100	
Shipping Weight	29 lbs.	29 lbs.	

Options: 75 different vinyl colors - see upholstery catalog

Foot Ring as an Accessory: Compatible on all 5 leg stools

Part#	Description	Shipping Weight
7901243000	Foot Ring for Leg Joint Comfort	6 lbs.



Seats: 15" diameter

Construction Specifications

Dimensions Polyfoam padding: 3" thick Height Range: 18" - 33" Materials **Bases:** Aluminum or black composite **Seats:** 75 colors of upholstered vinyl covering polyfoam padding. Clifton/Passaic Stool: 304 stainless steel Clifton Stool Legs: 18GA; Clifton Stool Seat: 14 GA Gauge Passaic Stool Legs: 18 GA Clifton/Passaic Stool: All welding by heliarc method and exposed welds Welding are ground smooth Finish Clifton/Passaic Stool: Exposed surfaces are free from pit marks, weld

Notes MR conditional stool available. See MR conditional catalog

For all vinyl colors, please see our upholstery catalog

2" hard rubber casters or 2" twin wheel plastic casters

seams, scale and are polished with #4 polish to a satin finish

Shipping Information

All shipping weights are estimates and may be changed without prior notice. Specifications may be improved and dimensions changed without prior notice. Please contact customer service for estimated shipping costs at info@blickman.com.

> Toll Free: (800) 247-5070 Local: (973) 330-0557

Standard

Casters



www.blickman.com info@blickman.com



MR Conditional

Exam Stools

EXAM STOOL

Description: Our MR conditional exam stool is industry standard, fully welded, stainless steel construction. Sit comfortably on 15" diameter, 3" thick padded seats. Choose between stationary feet tips or casters for mobility. Height is 21" on feet tips and 23" on MR casters.



Name	MR Conditional Stool	MR Conditional Stool with Casters
Item #	7445MR	7445MR-C
Model #	1027445000	1027445001
Shipping Weight	14 lbs.	16 lbs.

Construction Specifications

Standard Seats: 15" diameter **Dimensions** Polyfoam padding: 3" thick Height Range: 18" - 33" Materials Bases: Aluminum or black composite **Seats:** 75 colors of upholstered vinyl covering polyfoam padding. Clifton/Passaic Stool: 304 stainless steel Clifton Stool Legs: 18GA; Clifton Stool Seat: 14 GA Gauge Passaic Stool Legs: 18 GA Clifton/Passaic Stool: All welding by heliarc method and exposed welds Welding are ground smooth Finish Clifton/Passaic Stool: Exposed surfaces are free from pit marks, weld seams, scale and are polished with #4 polish to a satin finish Casters 2" hard rubber casters or 2" twin wheel plastic casters Notes MR conditional stool available. See MR conditional catalog For all vinyl colors, please see our upholstry catalog

Shipping Information

All shipping weights are estimates and may be changed without prior notice. Specifications may be improved and dimensions changed without prior notice. Please contact customer service for estimated shipping costs at info@blickman.com.



Foot Stools

Safe and sturdy to step up or down

Includes:

- 8" High Foot Stools
- 9" High Foot Stools
- Two Step Foot Stools
- EZ Stacking Stools + Dolly



THE BLICKMAN **DIFFERENCE**

- Whether chrome or stainless steel, our foot stools are high quality, durable, and safe.
- Every foot stool includes a non-slip rubber tread for secure footing and non-slip rubber tips on the stools' feet to secure floor placement.
- For heavy loads, some models include H-brace reinforcement.
- For additional user stability, some models include handrails or allow for the addition of a handrail.
- MR conditional options are also available.

Key Product Features



www.blickman.com



One Step

8" HIGH FOOT STOOLS

Description: Our MR conditional approved Kent foot stools include one step 8" off the ground. They are industry standard, fully welded, stainless steel construction with non-slip rubber treads and feet tips. 18" and 30" width sizes are available. Depth is 12". Include a 41" tall handrail for additional support.

8" stools have the same dimensions:

Width: 18" or 30" x Height: 8" x Depth: 12"



Name	18″W	30″W	
	With Handrail		
Model #	7757MR-HR	7761MR-HR	
Item #	1017757100	1017761100	
Shipping Weight	20 lbs.*	30 lbs.*	
	With	out Handrail	
Model #	7757MR	7761MR	
Item #	1017757000	1017761000	
Shipping Weight	11 lbs.	20 lbs.	

^{*}Shipping weight with an asterisk is the combined shipping weight of separate parts

Foot Stools

One Step

9" HIGH FOOT STOOLS

Description: Our 9" high foot stools come in either chrome or stainless steel and consist of one step 9" off the ground with non-slip rubber treads and feet tips. A 35 3/4" tall handrail for additional support is available for the chrome stool.

All 9" stools have the same dimensions:

Width: $14 \frac{1}{4}$ " x Height: 9" x Depth: $11 \frac{3}{8}$ "





Name	Stainless Steel Stool	Chrome Stool	
	With F	Handrail	
Model #	1261	1251	
Item #	1011261000	1011251000	
Shipping Weight	12 lbs.*	11 lbs.*	
	Without	t Handrail	
Model #	1260SS	1250	
Item #	1011260000	1021251640	
Shipping Weight	8 lbs.	7 lbs.	

^{*}Shipping weight with an asterisk is the combined shipping weight of separate parts



Two Step

TWO STEP FOOT STOOLS

Description: Our MR conditional approved two step foot stools are industry standard, fully welded, stainless steel construction with H-braces for additional reinforcement, non-slip rubber treads for secure footing, and feet tips for firm floor placement.



Step Surface Areas Top Step: 36" x 12" Bottom Step: 36" x 24" Heights of Steps Top Step: 15 3/4" Bottom Step: 8"

Name	2-Step Foot Stool	
	With Handrails (41" high)	_
Model #	7764MR-HR	
Item #	1012041096	
Shipping Weight	65 lbs.	

Two Step Foot Stool with Side Handrails



Step Surface Areas Top Step: 24" x 8" Bottom Step: 24" x 16" Heights of Steps Top Step: 15 3/4" Bottom Step: 8"

Name	The Donnelly
	With Handrails (41" high)
Model #	7763MR-HR
Item #	1017763100
Shipping Weight	40 lbs.*
	Without Handrail
Model #	7763MR
Item#	1017763000
Shipping Weight	27 lbs.

^{*}Shipping weight with an asterisk is the combined shipping weight of separate parts



EZ Stacking

OUR EZ Stacking Stool allows you to:

- ✓ Create wider platforms by locking side-to-side
- ✓ Create higher platforms by stacking top-to-bottom
- ✓ Store efficiently and transport easily with our corresponding dolly

EZ STACKING STOOLS



Description: Our EZ stacking stools and dolly are made with our durable stainless steel. The stools are covered with a non-slip rubber tread and feet tips.



EZ Stool Dimensions: 18.7/8" $\times 6.1/4$ " $\times 13.7/8$ "



Dolly Dimensions: $21 \frac{1}{2}$ " x 39" x 16 $\frac{1}{2}$ "

Name	EZ Stacking Stool	Dolly	
Model #	8862SS	8862SS-D	
Item #	1018862000	1018862300	
Ship Weight	14 lbs.	20 lbs.	



Construction Specifications

1 step: $6 \frac{1}{4}$ " - 9" high Standard **Dimensions** 2 step: Bottom step: 8" high, Top step: 15 3/4" high Materials Chrome 304 stainless steel Black non-slip rubber EZ Stacking Stools: 16GA platform, 14GA legs Gauge All Stainless Steel Stools: 16GA (Exception: 7757MR-HR and 7761MR-HR have 18GA handrails) Welding All welding by heliarc method and exposed welds are ground smooth Finish Exposed surfaces are free from pit marks, weld seams, scale and are polished with #4 polish to a satin finish Casters Dolly: 3" swivel casters with no brakes Notes MR conditional foot stools are available where indicated

Shipping Information

All shipping weights are estimates and may be changed without prior notice. Specifications may be improved and dimensions changed without prior notice. Please contact customer service for estimated shipping costs at info@blickman.com.

Note for: 1251 Stool

Parts shipped separately to avoid damage and lower costs. Shipping weight with asterisk is the combined shipping weight of separate parts. Some assembly required.



Hampers

For elimination of dirty linens

Includes:

- Round Hampers
- Square Hampers
- Space Saver Hamper
- Hamper Bags



THE BLICKMAN DIFFERENCE

- Our hampers facilitate easy and safe disposal of contaminated linens. They all move easily over floors on 2" or 3" casters. Our hampers have proven to be stable when empty or full due to their wide bases and added counterweights on covered models. A bottom brace holds and protects the hamper bag from contact with the floor.
- Our hampers come in either chrome or stainless steel and we offer different shapes and sizes to fit your needs and floor plans. Additionally, we offer covered models to control airborne pathogens. All covered hampers are foot-operated to help prevent cross contamination.
- For MR conditional environments, we offer two MR conditional models: one covered and one uncovered, for your choosing.
- We offer two optional hamper bags, nylon and water soluble, that fit our round hampers.


Key Product Features





ROUND HAMPERS

Description: We offer covered and uncovered round hampers in both stainless steel and chrome. They are 35" tall and roll on 2" hard rubber swivel casters. They come in two sizes 18" diameter and 25" diameter. The covered models operate via a spring-loaded lift mechanism that opens and closes the stainless steel lid. The pneumatic top allows for pin drop quiet operation.





Name Uncovered		Name	Cov	rered
Lid	Chrome	Stainless Steel	Chrome	Stainless Steel
Bag Size		18"	,	
Model #	2018	7774SS	2118	7774SS-LF
Item #	0962018000	0927774000	096211800P	092777420P
Shipping Weight	11 lbs.	11 lbs.	21 lbs.	25 lbs.
Bag Size		25″	•	
Model #	2025	8773SS	N/A	8773SS-LF
Item #	0962025000	0928773000	N/A	092877320P
Shipping Weight	12 lbs.	13 lbs.	N/A	35 lbs.

Hamper Bags: 18" and 25" white nylon hamper bags and water soluble hamper bags.



FOLDING ADJUSTABLE HAMPER

Description: Our folding chrome hamper adjusts to accommodate three hamper bag sizes with a bottom frame to support hamper bags. Travels on 2" hard rubber swivel casters.



Dimensions:

18" bag size: 19" x 36 1/2" x 18"
22" bag size: 19" x 40 1/2" x 22"
25" bag size: 19" x 42 1/4" x 25"

Name	Folding Adjustable Hamper
Model #	2230
Item #	0962230000
Shipping Weight	18 lbs.

SECURE CARE HAMPER

Description: Our secure care chrome hamper supports 36 - 42 gallon multipurpose bags. Its' polycoated steel lid opens using foot operated lift mechanism. Includes hamper storage pouch. Travels on 3" hard rubber swivels casters.



Dimensions:

• 18 5/8" x 37 3/4" x 19 1/2"

Name	Secure Care Hamper
Model #	2010
Item #	0962010000
Shipping Weight	26 lbs.





SPACE SAVER HAMPER

Description: Our space saver hamper offers hands free loading and unloading with a convenient non-slip foot pedal. Sturdy all stainless steel construction with bumper guards on lid and legs to protect equipment and walls. The solid bottom shelf provides extra security for heavy loads and a spring-loaded lift mechanism provides a pin-drop pneumatic top for quiet operation. Travels on 3" hard rubber swivel casters.



Name	Space Saver Hamper	
Model #	2301SS	
Item #	096230100P	
Shipping Weight	24 lbs.	



MR Conditional

Hamper

HAMPER

Description: Our MR conditional round stainless steel hamper travels easily on 2" MR casters and is compatible with our 18" nylon hamper bags and our water soluble hamper bags.





Name	MR Conditional Hamper	
Model #	7772MR	
Item #	0927772400	
Shipping Weight	20 lbs.	



Hampers

HAMPER BAGS



Disposable Water Soluble Plastic Hamper Bag



Reusable Nylon Hamper Bag



Secure Care Hamper Pouch

Part #	Description	Shipping Weight
98390000WS	100 Plastic Water Soluble Hamper Bags	9 lbs.
9839777201	Bag, 18" Diameter Ring Hamper Nylon White	1 lb.
9839777302	Bag, 25" Diameter Ring Hamper Nylon White	1 lb.
0962012000	Secure Care Hamper Pouch for Labels	1 lb.

Construction Specifications

Universal Dimensions	Hamper Bag Sizes: 18", 22", 25", 26"
Materials	304 stainless steel Chrome
Gauge	Stainless Steel: 1" diameter, 18GA tubing Chrome: 1" diameter tubing
Welding	All welding by heliarc method and exposed welds are ground smooth
Finish	Exposed surfaces are free from pit marks, weld seams, scale and are polished with #4 polish to a satin finish
Casters	2" & 3" hard rubber swivel casters (some with locks)
Notes	MR conditional hampers available where indicated

Shipping Information

All shipping weights are estimates and may be changed without prior notice. Specifications may be improved and dimensions changed without prior notice. Please contact customer service for estimated shipping costs at info@blickman.com.





Optimal continuous fluid delivery

Includes:

- IV Stands
- IV Stand Accessories
- IV Irrigation Safety Systems




THE BLICKMAN **DIFFERENCE**

IV Stands

- Our IV poles are optimal for maintaining continuous fluid delivery for every patient, available in a wide range of adjustable heights with a variety of easy locking options that can be operated with one hand.
- Our IV poles are all available in high quality chrome and stainless steel, and set on 4, 5, or 6 leg bases with casters, no projections, and great balance/tip resistance. Some bases are gravity weighted to help compensate for heavier loads. All bases have the ability to spin 360° in tight spaces.
- Most of our models have 2-, 4-, 6- and 8-(rake) hook adapters available. All come with colored secure grip tips.


Key Product Features





4 LEG CHROME IV STANDS

Description: We offer two 4-leg chrome IV stand models, one with a standard base, and one with a heavyweight base. Both base types include Wall Saver™ tips to protect against wall scuff marks. They both move on 2″ plastic twin wheel hooded swivel casters and they both come standard with two hooks. A four hook removable adapter is available for model 1305.



Name	Standard Base	Heavyweight Base
# of Hooks	2 Welded Hooks	2 Removable Hooks
Model #	1305	1310
Item #	0561305000	0561310000
Shipping Weight	13 lbs.	14 lbs.

4 Leg

4 LEG STAINLESS STEEL IV STANDS

Description: We offer two types of 4-leg stainless steel IV stands, one with a standard base, and one with a heavyweight base. They come standard with either 2 or 4 removable hooks. Contact customer service for 6- or 8-(rake) hooks.



Name	Standard Base	
# of Hooks	2 Hooks	
Model #	1410SS	
Item #	0571410000	
Shipping Weight	13 lbs.	
# of Hooks	4 Hooks	
Model #	1410SS-4	
Item #	0571410400	
Shipping Weight	14 lbs.	

Height Adjustment: 51 3/4" - 94"

Mechanism: Tru-Loc™

Base Material: Stainless steel with Wall Savers™, 22" diameter

Casters: 2" hard rubber swivel casters







Name Heavyweight Base

Lock Type	Twist Lock	Thumb Control
# of Hooks		2 Hooks
Model #	8890SS	8889SS
Item #	0528890000	0518889000
Shipping Weight	24 lbs.*	24 lbs.*
# of Hooks		4 Hooks
Model #	8890SS-4	8889SS-4
Item #	0528890400	0518889400
Shipping Weight	24 lbs.*	24 lbs.*

Height Adjustment: 66" - 100"

Mechanism: Twist Lock, Thumb Control, Foot Operated

Base Material: Heavyweight stainless steel flat bar, 22" diameter

Casters: 2" roller blade swivel casters

^{*}Shipping weight with an asterisk is the combined shipping weight of separate parts



5 Leg

5 LEG IV STANDS

Description: Our 5-leg IV stands come in either chrome (standard & heavy duty) or stainless steel poles and stand on a 24" diameter, polished die cast aluminum base with Wall Saver™ caps. Contact customer service for 6- or 8-(rake) hooks on certain models.

	Name	Chro	me Pole
	Pole	Standard (1")	Heavy Duty (11/4")
01/20	# of Hooks	2 Hooks	4 Hooks
	Caster Heights	2 1/2" Plastic	3" Rubber
	Height Adjustments	52 '1/ ₂ " - 93 ¹ / ₂ "	56" - 100"
100	Model #	1315-2	1350-4
1	Item #	0561315000	0561350400
	Shipping Weight	16 lbs.*	17 lbs.*
	# of Hooks	4 Hooks	6 Hooks
	Caster Heights	2 1/2" Plastic	3" Rubber
	Height Adjustments	52 '1/ ₂ " - 93 ¹ / ₂ "	56" - 100"
	Model #	1315-4	1350-6
	Item #	0561315400	0561350600
2 2	Shipping Weight	16 lbs.*	17 lbs.*

Height Adjustment Mechanism: Tru-Loc™ Base Material: 5 Leg, Polished die cast aluminum, 24" diameter Casters: 3" hard rubber OR 2 ½" twin wheel swivel casters








Name

Stainless Steel Pole



Height Adjustment Mechanism: Tru-Loc™, Twist Lock, Thumb Control Base Material: 5 Leg, Polished die cast aluminum, 24" diameter Casters: 2" or 3" hard rubber swivel casters

*Shipping weight with an asterisk is the combined shipping weight of separate parts

Name

6 LEG IV STANDS

Description: Our 6-leg IV stands are designed for hospital applications requiring heavier loads by using low center of gravity bases. They come in either a chrome or a stainless steel pole with a Tru-Loc™ height adjustment mechanism that allows the pole to reach a minimum of 52" high off the ground and a maximum of 90" off the ground. Our 6 leg IV stands come standard with 4-, 6-, and 8-(rake) removable hooks.



Features Include:

- √ 23 3/8" diameter base with low center of gravity for improved stability
- ✓ Heavy duty 3" roller blade swivel casters for easy maneuverability
- Steel powder coated finished base for durability (many colors available)

SS Pole

Accommodates up to 4 infusion pumps and up to 8 IV solution bags while preventing tubing line entanglement





Powder Coated Base


Stainless Steel Base

SS Pole

1371SS-4

0541371400

26 lbs.*



Chrome Pole



# of Hooks	6 Hooks		
Model #	1370-6	1370SS-6	137ISS-6
Item #	0561370600	0541370600	0541371600
Shipping Weight	29 lbs.*	29 lbs.*	26 lbs.*
# of Hooks		8 (rake) Hooks	
Model #	1370-RH	1370SS-RH	1371SS-RH
Item#	0561370RH1	05413700RH	05413710RH
Shipping Weight	29 lbs.*	29 lbs.*	26 lbs.*


^{*}Shipping weight with an asterisk is the combined shipping weight of separate parts

Custom Powder Coated Colors Available Upon Request

Please email us at info@blickman.com

5 LEG WASHABLE IV STAND

Description: Our Hackensack 5-leg IV stand is cart washable - built specifically to withstand high temperatures and cleaning fluids with weep holes in the pole and in the base for fluid drainage. This model comes with a stainless steel pole with a twist lock height adjustment mechanism (74" to 110" inches). It has a 22 7/8" diameter black composite base with twin wheel hooded swivel casters.




MR CONDITIONAL IV STAND

Description: Our 5 leg MR conditional IV stand consists of a stainless steel pole on a 24" diameter, polished diecast aluminum base with Wall Saver™ caps and 2" MR casters. It is height adjustable from 74" to 110" via a twist lock adjustment mechanism. Choose from either 2- or 4-hooks or contact customer service for 6 or 8 (rake) hooks.

# of Hooks	2 Hooks
Model #	7792MR-2
Item #	0537792002
Shipping Weight	16 lbs.*
# of Hooks	4 Hooks
Model #	7792MR-4
Item #	0537792004
Shipping Weight	17 lbs.*

MR Push Handle for MR IV Stands





Options

IV Hook Configurations









Part #	Description	Shipping Weight
780063700L	2 hooks	2 lbs.
780063800L	4 hooks	2 lbs.
780066900L	6 hooks	2 lbs.
7802808750	8 (rake) hooks	6 lbs.

Lock Options








Thumb Lock	Tru-Lock™	Twist Lock
------------	-----------	------------

Base Options

4 Leg Chrome Bases









Standard

Heavyweight

5 Leg Bases








Powder Coated

Stainless Steel



Accessories

IV POLE ATTACHMENTS



Accessory Clamp

Drainage Bag Hook Pump Support Tray Multi-pump Support

Part # 3501389001 3501381000

9101107600 7801380000

Description	Shipping Weight
Wire Basket	9 lbs.
Guidance Handle (MR Option Available)	4 lbs.
Power Strip	4 lbs.

Drainage Bag Hook	1 lbs.
Pump Support Tray	2 lbs.
Multi-pump Support	2 lbs.
Tank Holder	5 lbs.

Tank holder

4 lbs.

Construction Specifications

Universal Dimensions	Height Range: 47" - 110" Base Diameter Range: 22" - 24"
Materials	304 stainless steel Chrome
Gauge	Stainless Steel IV Stands: Base Pole: 1" Diameter, 16GA; Extension Pole: ¾" Diameter, 18GA Exception: 1370/1371 series Base Pole: 1 ¼" Diameter, 16 GA; Extension Pole: 1" Diameter, 18GA
Welding	All welding by heliarc method and exposed welds are ground smooth
Finish	Exposed surfaces are free from pit marks, weld seams, scale and are polished with #4 polish to a satin finish
Casters	See page 9
Notes	MR conditional IV stand option available where indicated

Shipping Information

All shipping weights are estimates and may be changed without prior notice. Specifications may be improved and dimensions changed without prior notice. Please contact customer service for estimated shipping costs at info@blickman.com.

Note for: IV Stands


Parts shipped separately to avoid damage and lower costs. Shipping weight with an asterisk is the combined shipping weight of separate parts. Some assembly required.

THE BLICKMAN **DIFFERENCE**

IV Irrigation Safety Systems

- Our IV Irrigation Safety Systems (IVTs) offer superior safety, efficiency, and financial value for fluid intensive IV irrigation procedures, such as arthroscopic, hysteroscopic, endoscopic, TURP, etc.
- They maintain continuous fluid delivery by offering infinitely adjustable heights for up to eight 3 liter or 5 liter IV bags on four independent channels.
- Their easy, one handed operation allows nurses to attach IV bags at comfortable heights, starting at 5 ft., and then raise bags up to 8 ft. (or to 9 ft. with custom extender) simply by pushing down on "squeeze/release" handles, an ergonomic solution that meets OSHA standards and greatly reduces the threat of back injuries.

Key Product Features



WHY CHOOSE OUR IVT?



Safe:

- This ergonomic solution meets OSHA standards and greatly reduces the threat of back injuries with its one-handed operation which allows nurses to attach IV bags at comfortable heights, simply by pushing down on "squeeze/release" handles, an
- Hard anodized aluminum construction with stain resistant finishes makes it so that no surfaces can chip, flake, or rust to introduce particle matter into a sterile or surgical environment.
- The lowest positions are cushioned by shock absorbers, thus sparing the common banging known to occur with the standard IV poles, which results in tearing of bag tabs and spilling remaining fluids onto the O.R. floor.
- The steady, five-legged bases resist potential tipping hazards, eliminating the need for "sand bagging".

Efficient:

- There are no electrical aspects that require overnight charging, so they are always ready, whenever you are ready.
- The ability to stagger heights of up to eight 3 liter or 5 liter bags ensures continuous flow of liquids, as it allows you to change an IV bag, while the other bags are still irrigating.
- Easy to clean and maintain
- Compact footprint frees up valuable space on the O.R. floor.

Affordable:

 Our IV Irrigation systems are value leaders among a few comparable models, saving you up to \$50/procedure on IV tubes vs. those required with far more expensive, pump driven systems.



IVTs

Description: We offer two types of IVTs, one with 2 hooks, which holds up to 4 IV bags and one with 4 hooks, which holds up to 8 IV bags. Hooks can raise and lower from 67" to 81". They travel on black plastic twin wheel hooded swivel casters. Both IVTs come with 26 1/2" diameter inverted, 5 leg base.

Our IVT stands have the same dimensions:

Width: 12 ³/₄" x Height: 89" x Depth: 12 ³/₄" Note: these dimensions do not include the base



Name	IVT 2 Channels	IVT 4 Channels
Item #	IVT-2	IVT-4
Model #	0541110200	0541110100
Shipping Weight	85 lbs.*	85 lbs.*

For Custom Extender

Please email us at info@blickman.com

*Shipping weight with an asterisk is the combined shipping weight of separate parts

Irrigation Systems

Construction Specifications

Construction Specifications

Universal Tower: Width: $12 \frac{3}{4}$ " x Height: 89" x Depth: $12 \frac{3}{4}$ "

Dimensions Base Diameter: $26 \frac{1}{2}$

Height Range of Hooks: 67" - 81"

Materials Hard anodized aluminum

Gauge N/A

Welding N/A

Finish Stain resistant finishes

Casters 2" black plastic, twin wheel, hooded swivel casters

Notes

Shipping Information

All shipping weights are estimates and may be changed without prior notice. Specifications may be improved and dimensions changed without prior notice. Please contact customer service for estimated shipping costs at info@blickman.com.

Note for: Parts shipped separately to avoid damage and lower costs.

IV Stands Some assembly required.



Mayo Stands

Maximum utility in minimal space

Includes:

- Chrome Mayo stands
- Stainless steel Mayo stands



THE BLICKMAN **DIFFERENCE**

- Our Mayo stands have evolved to meet the needs of the marketplace for maximum utility in an increasingly crowded surgical environment.
- We have an array of options with varying material, bases, shapes, and adjustment mechanisms to suit your needs.
- All of our Mayo stands are mobile and use 2" roller blade, hard rubber, or twin wheel plastic hooded casters.

Key Product Features





Mayo Stands

CHROME MAYO STANDS

Description: Each of our chrome Mayo stands contain a 12 /8" x 19 1/8" stainless steel tray and are height adjustable using the Tru-Loc™ friction knob. We offer different base and pole set ups to meet your varying needs.

Middle Pole

Base: Five leg cast aluminum base with wall savers™ and 2" twin wheel hooded casters

Double Pole

Base: Chrome loop base with two 2" twin wheel hooded back casters

Single Pole

Base: Two pronged chrome base with single hard rubber back casters







Name	Middle Pole	Double Pole	Single Pole
Height Adjustment	29 ¹ / ₂ " - 49"	34 ¹ / ₂ " - 54"	34 ¹ / ₂ " - 54"
Model #	1501	1510	1515
Item #	0661501000	0661510000	0661515000
Shipping Weight	14 lbs.	18 lbs.	12 lbs.



NEWARK MAYO STANDS

Description: Each of our stainless steel Mayo stands contain a 12 /8" x 19 1/8" stainless steel tray and are height adjustable using the Tru-Loc™ friction knob. We offer different base set ups to meet your varying needs.

Lightweight

Standard

Four Caster

Base: Two pronged stainless steel base with two hard rubber back casters

Base: Two pronged stainless steel base with two roller blade back casters

Base: Two pronged stainless steel base with four roller blade casters







Name	Lightweight	MR Stan	dard 🗽	Four Caster
Height Adjustment	37 1/4" - 64"	39 1/2" - 62"	391/2" - 62"	39 1/2" - 62"
Model #	1520SS	8841SS	8841MR	8848SS
Item #	0671520000	0638841000	0638841200	0648848000
Shipping Weight	26 lbs.	36 lbs.	36 lbs.	39 lbs.



BENJAMIN MAYO STAND

Description: Our Benjamin Mayo stand design is made for ease and durability. Our state of the art double post design minimizes lateral movement and increases overall stability. We offer this design in three different tray sizes.



Name



Via a foot operated mechanism, you can adjust the height from 39 1/2" to 62" high off the ground. A guard over the foot prevents accidental release.

Nume	benjamin mayo stanas	
Tray Size	12 ⁵ / ₈ " x 19 ¹ / ₈ " (Standard)	
Model #	8867SS	
Item#	0668867000	
Shipping Weight	40 lbs.	
Tray Size	16" x 21"	-
Model #	8869SS	
Item#	0668869000	
Shipping Weight	41 lbs.	
Tray Size	20" x 25"	
Model #	887ISS	
Item#	0668871000	
Shipping Weight	43 lbs.	

Reniamin Mayo Stands

Construction Specifications

Universal Dimensions	Standard Tray Size: Width: 12 5/8" x Depth: 19 1/8"
Materials	304 stainless steel Chrome
Gauge	All contain stainless steel trays: 20GA Benjamin Mayo Stands: ¾" diameter, 18GA tubes, 1/8" base Newark Mayo Stands: 1520SS: 1 ¼" diameter, 18GA tube, weighted base 8841SS/8848SS: 1" diameter, 16GA tube, 3/8" base
Welding	All welding by heliarc method and exposed welds are ground smooth
Finish	Exposed surfaces are free from pit marks, weld seams, scale and are polished with #4 polish to a satin finish
Casters	2" plastic hooded twin wheel, hard rubber, or roller blade casters
Notes	MR conditional Mayo stands available where indicated

Shipping Information

All shipping weights are estimates and may be changed without prior notice. Specifications may be improved and dimensions changed without prior notice. Please contact customer service for estimated shipping costs at info@blickman.com.



Hands free multi-station operation

Includes:

- Windsor Scrub Sinks
- Lodi Scrub Sinks
- Accessories




THE BLICKMAN **DIFFERENCE**

Scrub Sinks

- Our two scrub sink models, the Windsor and the Lodi, are each available with single, dual, and triple stations. A four-station unit is also available in the Windsor model.
- They have deep sink bowls with sound-deadened enclosed double walled panels that are angled in to divert water splashes along with spouts positioned toward the rear of the sink. Standard gooseneck faucets included.
- All operation options are hands free. The Lodi scrub sinks have infrared water control and knee-action water and soap dispensers, and the Windsor scrub sinks come standard with knee-action control, but they are also available in foot, elbow, and infrared sensor controls.

Key Product Features



Windsor

WINDSOR SCRUB SINKS

Description: Our Windsor scrub sinks are all wall mounted and utilize one undivided sink bowl. Choose from our four different hands free operations: knee, foot, elbow, or infrared. Standard gooseneck faucets included.

Inside Sink Bowl Dimensions

Width: 21 $\frac{3}{4}$ " x Height: 8" x Depth: 14 $\frac{1}{2}$ "



Name	1 Station	2 Stations	3 Stations	4 Stations
		Knee Action		
Model #	7878SS-KA	7879SS-KA	7880SS-KA	7881SS-KA
Item #	1317878000	1317879000	1317880000	1317881000
Shipping Weight	65 lbs.	145 lbs.	190 lbs.	305 lbs.
		Foot Action		
Model #	7878SS-FO	7879SS-FO	7880SS-FO	7881SS-FO
Item #	1317878003	1317879003	1317880003	1317881003
Shipping Weight	65 lbs.	145 lbs.	190 lbs.	305 lbs.
		Elbow Action		
Model #	7878SS-EA	7879SS-EA	7880SS-EA	7881SS-EA
Item #	1317878004	1317879004	1317880004	1317881004
Shipping Weight	65 lbs.	145 lbs.	190 lbs.	305 lbs.
		Infrared		
Model #	7878SS-IF	7879SS-IF	7880SS-IF	7881SS-IF
Item #	1317878001	1317879001	1317880001	1317881001
Shipping Weight	65 lbs.	145 lbs.	190 lbs.	305 lbs.

Lodi

LODI SCRUB SINKS

Description: Our Lodi scrub sinks come standard with both infrared water control and knee-action water and soap dispensers. All units are available in either pedestal or wall-mounted configurations. Standard gooseneck faucets included.



Width 32"



Width 64"



Width 96"



Height 53"

Name	1 Station	2 Stations	3 Stations
		Pedestal	
Model #	9881SSP	9882SSP	9883SSP
Item #	1339881P00	1339882P00	1339883P00
Shipping Weight	275 lbs.	420 lbs.	600 lbs.



Height 39 1/4"

Wall Mounted				
Model #	9881SSW	9882SSW	9883SSW	
Item #	1339881W00	1339882W00	1339883W00	
Shipping Weight	215 lbs.	360 lbs.	570 lbs.	

Accessories

LODI SCRUB SINK ACCESSORIES

- Factory Installed Options -





Digital Timer



Eyewash

- For Wall Mount Installation -



Chair CarrierFor Wall Mount Installation

of Stations One Two Three

Part # 4301224751 4301224752 4301224753

Shipping Weight 47 lbs. 119 lbs. 141 lbs.



Construction Specifications

Construction Specifications

Standard **Lodi:** (W) 32", 64", 96" x (D) 27 1/4 x (H): 39 1/4", 53"

Dimensions Internal Sink (Single Station): 28" x 18 1/2" x 20 1/4"

Windsor: (W) 25", 47", 69", 91 x (D) 18 1/2" x (H): 15 3/4"

Internal Sink: 21 3/4" x 8" x 14 1/2"

Materials/ 304 stainless steel

Gauge Lodi: Mix of 16GA and 14GA

Windsor: 16 GA

Double/Triple Station Lodi Scrub Sinks: Plastic Dividers

Welding All welding by heliarc method and exposed welds are ground smooth

Finish Exposed surfaces are free from pit marks, weld seams, scale and are

polished with #4 polish to a satin finish

Electrical Available in standard 100/120 VAC input or 220/240 VAC

Units are FDA registered and conform to standards for Safety Medical Electrical Equipment UL 606011 (US) and CSA C22.2 #601.1M90 (Canada)

Units are cord connected with a hospital grade cap

Casters N/A

Shipping Information

All shipping weights are estimates and may be changed without prior notice. Specifications may be improved and dimensions changed without prior notice. Please contact customer service for estimated shipping costs at info@blickman.com.





Solution Stands/ Receptacles

Roll, collect, dispose

Includes:

- Solution Stands
- Kickbucket
- Receptacle



THE BLICKMAN **DIFFERENCE**

- Our high quality, heavy-gauge, noise-reducing, solution stands and receptacles ensure your solutions and waste products are easily collected and disposed of in these easy-to-move, easy-to-clean product.
- Our solution stands and receptacles move effortlessly and quietly on 2" swivel casters, and are made out of our high-quality, easy-to-clean 304 stainless steel with MR conditional options available for both.
- Our solution stands are 33" high for compatibility with surgical table heights and are in the safe zone for reduced back stress, and our receptacles are low to the floor to pass under surgical tables and make it easy to drop waste into without reaching around the table.
- Our basins and pails are removable for both the solution stands and receptacles that are held in place by rubber bumpers that also help prevent noise from metal on metal contact.

Key Product Features



Solution Stands

Description: We offer two types of stainless steel solution stands: a single basin model and a double basin model. Our solution stands are 33" tall offering compatibility with surgical table heights and are in the safe zone for reduced back stress. Choose between a lower shelf or an H-brace for added storage and support. MR conditional models available.





Dimensions:Width: 15" x Height: 33" x Depth: 15"

Dimensions:Width: 29" x Height: 33" x Depth: 15"

		Name	Single Basin Solution Stand		Double Basin S	Solution Stand
	-		(MR)	MR Sh	elf MR	MR
		Model #	7807SS	7807MR	7808SS	7808MR
Į.	-W	Item #	0717807000	0717807066	0727808000	0727808066
6	6	Shipping Weight	20 lbs.	20 lbs.	33 lbs.	33 lbs.
				Н-В	race	
	- `	Model #	7807SS-HB	7807MR-HB	7808SS-HB	7808MR-HB
	-	Item #	0717807100	0717807166	0727808100	0727808166
P	P	Shipping Weight	19 lbs.	19 lbs.	24 lbs.	24 lbs.

Solution Basin Also Sold Separately

Name	Model #	Item #	Ship Weight
8.5 qt. Solution Basin Only	SB-1	9737370120	3 lbs.

Receptacles

Description: Our stainless steel kickbucket consists of a continuous encircling rubber bumper, which helps prevent tipping, wall scuff marks, and noise. A stainless steel cross brace provides excellent protection, securely holds a heavy gauge 13 qt. pail, as well as roller blade casters for mobility.

Parts Also Sold Separately





#KB-1



#7766NP

Parts Also Sold Separately





#SB-1



#8766NP

Kickbucket

13 qt. Pail

Receptacle

8.5 qt. Solution Basin

KICKBUCKET

Name	Model #	Item#	Ship Weight
Lenox Kickbucket with Pail	7766SS	0817766000	12 lbs.
Lenox Kickbucket with Pail	7766MR	0817766066	12 lbs.
	ORDER SEPARAT	ELY	
Lenox Kickbucket Frame Only	7766NP	0817766001	9 lbs.
13 qt. Pail Only	KB-1	9701920000	4 lbs.

RECEPTACLE

Name	Model #	Item#	Ship Weight
Sponge Receptacle with Basin	8766SS	0828766000	11 lbs.
	ORDER SEPARAT	ELY	
Sponge Receptacle Frame Only	8766NP	0828766001	9 lbs.
8.5 qt. Solution Basin Only	SB-1	9737370120	3 lbs.



Construction Specifications

Universal **Solution Stands:** 15", 29" x 33" x 15"

Dimensions Solution Basin: 8.5 qts.

Pail: 13 qts.

Materials 304 stainless steel

Gauge Solution Stand Frame: 1" diameter tube, 16GA

Solution Stand Shelf: 20GA Pail: 22GA; Basin: 14GA

Welding All welding by heliarc method and exposed welds are ground smooth

Finish Exposed surfaces are free from pit marks, weld seams, scale and are

polished with #4 polish to a satin finish

Casters 2" swivel casters

Notes MR conditional options available where indicated

Shipping Information

All shipping weights are estimates and may be changed without prior notice. Specifications may be improved and dimensions changed without prior notice. Please contact customer service for estimated shipping costs at info@blickman.com.





Safe, sterile, stainless steel

Includes:

- **Dressing Jars**
- Solution Basins Instrument Trays
- Emesis Basins •
- Ointment Jars **Bucket Pail**
- - **lodine Cups**

Catheter Tray

- Medicine Cup
- Graduated Measuring Cups



THE BLICKMAN **DIFFERENCE**

- We offer a comprehensive range of surgical utensils to complement the high-end quality of our hospital equipment.
- We offer utensils that function to hold other surgical materials and liquids of various types.
- All utensils are stainless steel and welded smooth to ensure sterility for surgery.
- Our measuring utensils have units of measurement commonly used in the medical industry etched inside.

Key Product Features



Description: Our surgical utensils are solid stainless steel and are great for sterile use in the OR.



Model #	Item#	Capacity/Dimensions	Shipping Weight			
		Dressing Jars				
DJ-1	9788000001	0.5 liters	1 lbs.			
DJ-2	9788010002	1.0 liters	1 lbs.			
DJ-3	9788020003	1.9 liters	2 lbs.			
DJ-4	9788030004	2.8 liters	2 lbs.			
DJ-5	9788040005	4.0 liters	2 lbs.			
DJ-6	9788060006	5.7 liters	2 lbs.			
OJ-7	9788080007	7.6 liters	3 lbs.			
		Emesis Basins				
EB-1	9717090400	12 oz.	1 lbs.			
EB-2	9720595450	16 oz.	1 lbs.			
EB-3	9725511555	28 oz.	1 lbs.			
EB-4	9730514065	56 oz.	56 oz. 1 lbs.			
		Ointment Jar				
OJ-1	9787720001	8 oz.	1 lbs.			



Bucket Pail # KB-1



Solution Basins

Model #	Item#	Capacity/Dimensions	Shipping Weight
		Bucket Pail	
KB-1	9701920000	12.3 liters	4 lbs.
		Solution Basin	
SB-1	9737370120	8 liters	3 lbs.
1830	9701830000	6.5 liters	2 lbs.

Description: Our surgical utensils are solid stainless steel and are great for sterile use in the OR.





Model#	Item#	Capacity/Dimensions	Shipping Weight
		Instrument Trays	
IT1015	9710015001	15" W x 10" D	3 lbs.
IT1621	9716002100	21" W x 16" D	4 lbs.
IT1525	3603915250	25" W x 16" D	5 lbs.
IT2025	3603920251	25" W x 20" D	6 lbs.
IT1219	9719126250	19" W x 12 ¹ / ₂ " D	4 lbs.
IT2520	9736392025	20" W x 25" D	6 lbs.
		Catheter Tray	
CT-1	9785030101	8 ½" W x 3" D x 1½" H	1 lbs.
CT-IT	9782831000	8 1/₂" W x 3" D x 1" H	1 lbs.
CT-1C	9782832000	8 ¹ / ₂ " W x 3" D x 1 ¹ / ₂ " H	1 lbs.





Medicine Cup
MC-1



Graduated
Measuring Cups
GM1-3

Model #	Item #	Capacity/Dimensions	Shipping Weight		
		lodine Cups			
IC-1	9775500002 6 oz.		1 lbs.		
IC-2	9710565001	14 oz.	1 lbs.		
		Medicine Cup			
MC-1	9784920001	2 oz.	1 lbs.		
	Gra	duated Measuring Cups			
GM-1	9795160001	16 oz.	1 lbs.		
GM-2	9795320002	32 oz.	1 lbs.		
GM-3	9795640003	64 oz.	1 lbs.		

Construction Specifications

	Construction Specifications
Standard Dimensions	N/A
Material	304 stainless steel
Gauge	 20GA: Graduated Measuring Cup, Instrument trays (exception 15"x10" Instrument Tray is 24GA) 22GA: Dressing Jars, Emesis Basins, Ointment Jar, Bucket Pail, Solution Basin, Catheter Tray and Cover, Iodine Cups, Graduated Measuring Cups
Welding	All welding by heliarc method and exposed welds are ground smooth
Finish	Exposed surfaces are free from pit marks, weld seams, scale and are polished with #4 polish to a satin finish
Casters	N/A

Shipping Information

All shipping weights are estimates and may be changed without prior notice. Specifications may be improved and dimensions changed without prior notice. Please contact customer service for estimated shipping costs at info@blickman.com.



Adaptable work and prep surfaces

Includes:

- Instrument Tables
- Utility Tables
- Work Tables




THE BLICKMAN **DIFFERENCE**

Adjustable Tables

- Our adjustable tables are instrument tables designed to straddle the OR table using either manual or motorized adjustable height mechanisms to accommodate changes in the table position during neuro, cardiac, and orthopedic surgery, as well as varying heights of personnel during general surgery.
- You can customize the tabletop dimensions to fit your differing surface area needs and the leg clearance dimensions to fit your differing OR table dimensions.
- Models are made out of stainless steel and include casters for mobility.
- The motorized table moves along telescoping aluminum lifting columns.

Key Product Features






BACK TABLE

Description: Our specialty back table is designed to provide additional room for extensive surgeries which require a wide array of equipment and instruments. Height adjustment from 46-49" is achievable with either pneumatic assist or manually. The top tier 'hide-away' shelf is angled for easier instrument viewing and holds multiple trays. Tables have sanitary closed tops which help prevent bacterial build up. Sealed reinforcing channels provide rigidity and durability.

Performance Back Table





Manual Back Table



Model#	Item#	Description	Length	Depth	Height	Overall Adjustable Height	Shipping Weight
BT4830	0364830340	Back Table with 3" Overshelf	48"	30"	34"	46-49"	250 lbs.
BT6030	0366030340	Back Table with 3" Overshelf	60"	30"	34"	46-49"	298 lbs.
BTP6030	036603034P	Back Table Performance with 3" Overshelf	60"	30"	34"	46-49"	360 lbs.
BT7230	0367230340	Back Table with 3" Overshelf	72"	30"	34"	46-49"	345 lbs.
BTP7230	036723034P	Back Table Performance with 3" Overshelf	72"	30"	34"	46-49"	413 lbs.

^{*}All models come standard with 3" casters.

Custom Dimensions Available Upon Request


Please email us at info@blickman.com



THE LIBERTY TABLE

Description: Our Liberty table was designed to straddle the OR table with adjustable height mechanisms. It travels effortlessly on antistatic swivel casters with brakes and utilizes a motorized height adjustment mechanism. Use a foot pedal or hand remote to operate the motorized adjustment. Two standard surface areas available.

Customize the surface area dimensions to fit your needs.



Standard Features





Foot Pedal

Hand Remote

MOTORIZED ADJUSTMENT:

Use the motorized foot operated/hand operated height adjustment mechanism.

Model #	Item#	Description	Width	Depth	Clearance	Shipping Weight
7902SS-M	0157902000	Liberty Table with Adjustable Height 38 3/4-58"	41"	46"	29"	315 lbs.
7901SS-M	0157901000	Liberty Table with Adjustable Height 38 3/4-58"	50 1/2	" 24"	38 3/8"	290 lbs.

Custom Dimensions Available Upon Request

Please email us at info@blickman.com



HYDRAULIC INSTRUMENT TABLE

Description: Our hydraulic instrument table is height adjustable to ensure compatibility with multiple OR table heights and healthcare personnel heights. Adjust its height with a manual hand crank from 30" - 42". It includes one lower shelf and stands stationary on glider feet tips. Select the surface area with dimensions that fit your needs.

Customize the surface area dimensions to fit your needs.



Model#	Item#	Description	Length	Width	Shipping Weight
7835HYD	0177835000	Hydraulic Table with Adjustable Height 30-42"	48"	24"	145 lbs.
7836HYD	0177836030	Hydraulic Table with Adjustable Height 30-42"	60"	30"	188 lbs.

Custom Dimensions Available Upon Request


Please email us at info@blickman.com



THE GEORGETOWN TABLE

Description: Our Georgetown table was designed to straddle the OR table with a manual hand crank adjustable height mechanism. Height adjustment within ergonomic range between 36" and 56" off the ground. Travels easily on 3" swivel casters with breaks. Two standard surface areas available. Select extended arm design for electrophysiology procedures.

Customize the surface area dimensions to fit your needs.



Name

7892SS-0

Georgetown Table

7895SS-0 CRT

Model#	Item#	Width	Depth	Leg Clearance	Shipping Weight
7892SS-0	0157892000	36"	20"	28 ⁵ / ₈ "	135 lbs.
7895SS-0	0157895000	48"	24"	40 5/8"	160 lbs.
7895SS-O CRT	0157895062	62"	24"	40"	190 lbs.
7896SS-0	0157896000	34"	46"	26"	175 lbs.

7895SS-0

Custom Dimensions Available Upon Request

Please email us at info@blickman.com

7896SS-0



OVERBED TABLE

Description: Our overbed table comes in an attractive mark resistant laminate top in a walnut wood grain. It moves up or down with a simple one touch adjustment. Security features include a built-in column lock and four point base. The table travels smoothly on four 2" hooded casters.

Overbed Table Dimensions

Width: 30" x Height: 30" (max) x Depth: 15 1/2"



Model#	Item#	Description	Width	Depth	Height	Shipping Weight	
3400	2833400000	Overbed Table	30"	15 1/2"	30" Max	27 lbs.	





Construction Specifications

Universal Dimensions	Adjustable heights to fit over OR table or hospital bed
Materials	304 stainless steel Liberty: Telescoping aluminum lifting columns. Overbed Table: Mar resistant laminate (tray), chrome plate steel (column & base)
Gauge	Liberty: 16GA Top; 11GA Sub Top; 14GA Base Georgetown: 18GA Top; 16GA Sub Top, 2" 16GA square tube legs
Welding	All welding by heliarc method and exposed welds are ground smooth
Finish	Exposed surfaces are free from pit marks, weld seams, scale and are polished with #4 polish to a satin finish
Casters	Liberty: 4" antistatic swivel casters with brakes Georgetown: 3" swivel casters, two with brakes
Notes	Liberty: Motorized adjustment mechanism Electrical Ratings: 100-120VAC, 15AMP, 60 Hz. Georgetown: Manual adjustment mechanism

Shipping Information

All shipping weights are estimates and may be changed without prior notice. Specifications may be improved and dimensions changed without prior notice. Please contact customer service for estimated shipping costs at info@blickman.com.

Custom Dimensions Available Upon Request

THE BLICKMAN **DIFFERENCE**

Stationary Tables

- We offer a variety of stationary tables to fit any healthcare need. Our stationary tables include instrument tables, utility tables, and general work tables. Some models also come in MR conditional equivalents.
- Our instrument tables are reinforced with H-braces for stability or include an optional lower shelf for storing supplies. We offer standard dimensions or that can be customized to any size. We also have an optional nested set of up to six tables that can be ordered together or separately.
- Our utility tables come in a variety of drawer and shelf configurations to store and protect
 equipment and can include a workspace so you can store and prep your supplies on one
 efficient and transportable, setup.
- Our work tables are heavy duty with large surface areas designed for use in sterile and non-sterile preparation.

Key Product Features





NESTED TABLES

Description: Our nested instrument tables are space saving alternatives with a unique design that requires a mere 8 square feet of floor space, but yields 31 square feet of usable work space. All tables travel on 3" casters.



Name	Dimensions (WxHxD)	Model #	Item#	Ship ping Weight
Complete Set of 6	8 Sq. Ft. of Floor Space 31 Sq. Ft. of Work Space	782ISS	0137821000	255 lbs.
	ORDE	R SEPARATELY		
Largest	48" x 42" x 24"	7822SS	0137822000	56 lbs.
Second	44" x 40" x 22"	7823SS	0137823000	52 lbs.
Third	40" x 38" x 20"	7824SS	0137824000	51 lbs.
Fourth	36" x 36" x 18"	7825SS	0137825000	37 lbs.
Fifth	32" x 34" x 16"	7826SS	0137826000	34 lbs.
Smallest	28" x 32" x 14"	7827SS	0137827000	32 lbs.



MR Conditional

Tables

NESTED TABLES

Description: Our MR conditional nested instrument tables are space saving alternatives with a unique design that requires a mere 8 square feet of floor space, but yields 31 square feet of usable work space. All nested tables travel on 2" MR casters.





Name	Dimensions (WxHxD)	Model #	Item#	Ship Weight
Complete Set of 6	8 Sq. Ft. of Floor Space 31 Sq. Ft. of Work Space	7821MR	0197821000	255 lbs.
	ORDER S	EPARATELY		
Largest	48" x 42" x 24"	7822MR	0197822100	56 lbs.
Second	44" x 40" x 22"	7823MR	0197823100	52 lbs.
Third	40" x 38" x 20"	7824MR	0197824100	51 lbs.
Fourth	36" x 36" x 18"	7825MR	0197825100	37 lbs.
Fifth	32" x 34" x 16"	7826MR	0197826100	34 lbs.
Smallest	28" x 32" x 14"	7827MR	0197827100	32 lbs.



HOWARD & SAWYER TABLES

Description: Our instrument tables offer complete 18 gauge, 304 stainless steel construction. They are lightweight, yet extremely durable. Sanitary closed top help prevents bacterial build up. Sealed reinforcing channels provide rigidity. All tables are 34" tall, which is in the safety zone for reduced back stress, and travel on 3" casters. There is a 16" clearance between shelves. Our 4 most common sizes are listed below, however we offer a wide variety of other sizes that you can order by calling our customer service or visiting our website.

Sawyer Table

Howard Table





Table	Sawyer Table		Howard Table			
Dimensions (WxD)		H-Brace			Shelf	
20" x 16"	7840SS	0127840000	28 lbs.	7830SS	0117830000	32 lbs.
30" x 16"	7839SS	0127839000	33 lbs.	7829SS	0117829000	40 lbs.
33" x 18"	7841SS	0127841000	39 lbs.	7831SS	0117831000	47 lbs.
36" x 20"	7842SS	0127842000	42 lbs.	7832SS	0117832000	52 lbs.
36" x 24"	7843SS	0127843000	48 lbs.	7833SS	0117833000	56 lbs.
48" x 20"	7844SS	0127844000	53 lbs.	7834SS	0117834000	66 lbs.
48" x 24"	7845SS	0127845000	56 lbs.	7835SS	0117835000	69 lbs.
60" x 24"	7846SS	0127846000	178 lbs.	7836SS	0117836000	190 lbs.
72" x 24"	7847SS	0127847000	186 lbs.	7837SS	0117837000	210 lbs.

Options: Custom dimensions, 40"H tables, three sided guardrails, and drawers (see Crescent table design).



MR Conditional

Tables

INSTRUMENT TABLES

Description: Our MR conditional instrument tables offer complete 18 gauge, 304 stainless steel construction. A sanitary closed top helps prevent bacterial build up. Sealed reinforcing channels provide increased rigidity. All tables are at an ergonomic reinforcing height of 34" tall with 1 ¼" diameter legs that are fully welded for strength and stability. Each MR conditional instrument table travels on 2" MR casters. There is a 16" clearance between shelves. Our 4 most common sizes are listed below, however we offer a wide variety of other sizes which you can order by calling customer service or visiting our website.



	_	_
Surface	Araa	

MR Conditional Sawyer Table

MR Conditional Howard Table

20" x 16" Model # 7840MR 7830MR Item # 0197840100 0187830100 36" x 24" Model # 7843MR 7833MR Item # 0197843100 0187833100 56 lbs. 56 lbs. 48" x 24" Model # 7845MR 7835MR Item # 0197845100 0187835100 Shipping Weight 56 lbs. 69 lbs. Model # 7846MR 7836MR Item # 0197846100 0187836100 Shipping Weight 189 lbs. 200 lbs.		•	
Item # 0197840100 0187830100 Shipping Weight 28 lbs. 32 lbs. 36" x 24" Model # 7843MR 7833MR Item # 0197843100 0187833100 56 lbs. 48" x 24" Model # 7845MR 7835MR Item # 0197845100 0187835100 Shipping Weight 56 lbs. 69 lbs. 60" x 24" Model # 7846MR 7836MR Item # 0197846100 0187836100			20" x 16"
Shipping Weight 28 lbs. 32 lbs. 36" x 24" Model # 7843MR 7833MR Item # 0197843100 0187833100 56 lbs. 48" x 24" Model # 7845MR 7835MR Item # 0197845100 0187835100 Shipping Weight 56 lbs. 69 lbs. 60" x 24" Model # 7846MR 7836MR Item # 0197846100 0187836100	Model #	7840MR	7830MR
36" x 24" Model # 7843MR 7833MR Item # 0197843100 0187833100 48" x 24" Model # 7845MR 7835MR Item # 0197845100 0187835100 Shipping Weight 56 lbs. 69 lbs. 60" x 24" Model # 7846MR 7836MR Item # 0197846100 0187836100	Item #	0197840100	0187830100
Model # 7843MR 7833MR Item # 0197843100 0187833100 Shipping Weight 48 lbs. 56 lbs. 48" x 24" Model # 7845MR 7835MR Item # 0197845100 0187835100 Shipping Weight 56 lbs. 69 lbs. 60" x 24" Model # 7846MR 7836MR Item # 0197846100 0187836100	Shipping Weight	28 lbs.	32 lbs.
Item # 0197843100 0187833100 Shipping Weight 48 lbs. 56 lbs. 48" x 24" Model # 7845MR 7835MR Item # 0197845100 0187835100 Shipping Weight 56 lbs. 69 lbs. Model # 7846MR 7836MR Item # 0197846100 0187836100			36" x 24"
Shipping Weight 48 lbs. 56 lbs. 48" x 24" Model # 7845MR 7835MR Item # 0197845100 0187835100 Shipping Weight 56 lbs. 69 lbs. 60" x 24" Model # 7846MR 7836MR Item # 0197846100 0187836100	Model #	7843MR	7833MR
48" x 24" Model # 7845MR 7835MR Item # 0197845100 0187835100 Shipping Weight 56 lbs. 69 lbs. 60" x 24" Model # 7846MR 7836MR Item # 0197846100 0187836100	Item #	0197843100	0187833100
Model # 7845MR 7835MR Item # 0197845100 0187835100 Shipping Weight 56 lbs. 69 lbs. 60" x 24" Model # 7846MR 7836MR Item # 0197846100 0187836100	Shipping Weight	48 lbs.	56 lbs.
Item # 0197845100 0187835100 Shipping Weight 56 lbs. 69 lbs. 60" x 24" Model # 7846MR 7836MR Item # 0197846100 0187836100			48" x 24"
Shipping Weight 56 lbs. 69 lbs. 60" x 24" Model # 7846MR 7836MR Item # 0197846100 0187836100	Model #	7845MR	7835MR
60" x 24" Model # 7846MR 7836MR Item # 0197846100 0187836100	Item #	0197845100	0187835100
Model # 7846MR 7836MR Item # 0197846100 0187836100	Shipping Weight	56 lbs.	69 lbs.
Item # 0197846100 0187836100			60" x 24"
	Model #	7846MR	7836MR
Shipping Weight 189 lbs. 200 lbs.	Item #	0197846100	0187836100
	Shipping Weight	189 lbs.	200 lbs.

Custom Dimensions Available Upon Request



ANGULAR INSTRUMENT TABLE

Description: Choose our angular instrument table to have separate left and right work spaces. Add a lower shelf or drawers for additional storage. This 34" tall table is in the safety zone for reduced back stress, and it travels on 3" casters. Select the surface area with dimensions that fit your needs.



Model#	Item#	Description	Width	Depth	Height	Shipping Weight	
7862SSA	0147862000	Angular Instrument Table	60"	15"	34"	120 lbs.	
7863SSA	0147863000	Angular Instrument Table	72"	18"	34"	137 lbs.	

Custom Dimensions Available Upon Request



UTILITY TABLES

Description: Our utility tables are popular in the OR as anesthesia tables and throughout healthcare facilities as a general utility table. Our superior construction adds a supporting channel underneath the top shelf to enhance its strength. All of our utility tables are 34" tall, which is in the safety zone for reduced back stress, and travel easily on 3" casters.

Shelf Configuration Includes:

- ✓ Two shelves, 17" and 34" off the ground
- ✓ Four sided guardrails for protection

Drawer Configurations Include:

- ✓ A positive stop to prevent accidental drawer removal
- ✓ Ball bearing slides for easy sliding drawers
- Three sided guardrails for protection (optional four sided guardrails are available)

All of our utility tables are the same dimensions:

Width: 20" x Height: 34" x Depth: 16"











Lower Shelf Clearance:

Lower Shelf Clearance:

Lower Shelf Clearance: Lower Shelf Clearance: Lowest Drawer Clearance: 4 1/4"

# of Drawers	0	1	2	3	4
Model #	7852SS	7850SS	7754SS	7857SS	7859SS
Item #	0257852000	0217850000	0227754000	0237857000	0247859000
Shipping Weight	32 lbs.	46 lbs.	59 lbs.	72 lbs.	82 lbs.

Options: Feet: glides, rubber tips, or larger diameter casters, four sided guard rails on drawer models



CRESCENT TABLES

Description: Our crescent table has the same construction as the Howard and Sawyer Instrument tables with the addition of drawers and guard rails. This table includes a large work surface area (36"x 20"), two 4" deep drawers and one optional lower shelf 16" off the ground. This model is 34" tall, which is in the safety zone for reduced back stress, and travels easily on 3" casters.

Customize the surface area dimensions of these models to fit your needs.





Model#	Item#	Description	Width	Depth	Height	Shipping Weight
7856SS	0317856000	Crescent Table (H-Brace)	36"	20"	34"	78 lbs.
7856SS-SH	0317856003	Crescent Table (Shelf)	36"	20"	34"	81 lbs.

Options: Feet: glides, rubber tips, or larger diameter casters



WORK TABLES

Description: Our work tables are designed for use as prep tables. They are heavy duty with large surface areas designed for use in sterile and non-sterile environments. The smallest table dimensions use four legs, and all other larger dimensions use six legs. All of our work tables are 36" tall, which is in the safety zone for reduced back stress., and come standard with bullet feet. Rubber tips, casters, and gliders are available as alternative options.

Customize the surface area dimensions of these models to fit your needs.



Model#	Item#	Description		Width	Depth	Height	Shipping Weight
			H-Brace				
CSDT6830	0336830000	Work Table		68"	30"	36″	210 lbs.
CSDT8036	0338036000	Work Table		80"	36"	36"	300 lbs.
CSDT9244	0339244000	Work Table		92"	44"	36"	390 lbs.
Shelf							
CSDT6830S1	0336830100	Work Table		68"	30"	36"	260 lbs.
CSDT8036S1	0338036100	Work Table		80"	36"	36"	355 lbs.
CSDT9244S1	0339244100	Work Table		92″	44"	36"	435 lbs.

Options: Feet: glides, rubber tips, or larger diameter casters

Custom Designs Available Upon Request



STATIONARY TABLE OPTIONS

Drawers

Description: For added storage, add drawers beneith any tabletop.



*contact us for available drawers and sizing

Guardrails

Description: For added placement security, add guardrails to any tabletop.



Custom Designs Available Upon Request



STATIONARY TABLE OPTIONS

Feet Options

Rubber Tips | Casters | Gliders









Custom Designs Available Upon Request



Construction Specification

Universal Dimensions	Custom dimensions available on all tables. Please email us at info@blickman.com.
Materials	304 stainless steel
Gauge	Instrument Tables: 111/4" tube legs, 18GA table top Utility Tables: 1" Diameter 18GA tube legs, 18GA top & body Work Tables: 15/8", 16GA tube legs, 16GA table top & sub top
Welding	All welding by heliarc method and exposed welds are ground smooth
Finish	Exposed surfaces are free from pit marks, weld seams, scale and are polished with #4 polish to a satin finish
Casters	3" Casters
Notes	MR conditional tables available where indicated

Custom Dimensions/Designs Available Upon Request

Please email us at info@blickman.com

Shipping Information

All shipping weights are estimates and may be changed without prior notice. Specifications may be improved and dimensions changed without prior notice. Please contact customer service for estimated shipping costs at info@blickman.com.

Toll Free: (800) 247-5070 Local: (973) 330-0557



www.blickman.com info@blickman.com



Utility Carts

Convenient transport for medical equipment

Includes:

- Light Duty
- Medium Duty
- Heavy Duty



THE BLICKMAN **DIFFERENCE**

- Our utility carts are ideal for transporting medical equipment which needs to be used in various areas in your medical facility. We offer light duty, medium duty, and heavy duty options to suit all of your needs.
- Common features include a push handle, open construction for easy access, guardrails for securing contents, and 4" swivel casters for easy maneuverability.
- They contain two or three solid stainless steel shelves to provide ample storage space.
- The stainless steel construction allows for easy cleaning and sanitizing, while resisting corrosion and stains.
- We also offer a folding utility cart with wire baskets for any application such as movement of mail, patient records, and laboratory materials.

Key Product Features



Utility Carts

Light Duty

FOLDING UTILITY CART

Description: Our folding utility cart is versatile enough for any application such as movement of mail, patient records, and laboratory materials. Choose any combination of 6" or 12" deep wire baskets to complete your cart. The metal frame is designed for easy assembly and folds for storage. Rubber stoppers support and protect the wire racks and provide a quiet ride on four, 3" hard rubber casters.

Folding Utility Cart Dimensions:

Width: 29 1/4" x Height: 38 1/2" x Depth: 19 1/2"





6" Wire Basket (2441)
Item #: 2422441000



12" Wire Basket (2442) Item #: 2422442000

Name	Model#	Item #	Shipping Weight				
Kit with 6" + 12" Baskets	2440 KIT	24224400KT	33 lbs.				
Kit with Two 6" Baskets	2440 6 KIT	24224401KT	31 lbs.				
Kit with Two 12" Baskets	2440 12 KIT	24224402KT	35 lbs.				
ORDER SEPARATELY							
Folding Frame	2440	2422440000	21 lbs.				
6" Wire Basket	2441	2422441000	5 lbs.				
12" Wire Basket	2442	2422442000	6 lbs.				

Utility Carts

Light Duty

LIGHT DUTY CARTS

Description: Our light duty 7534 utility cart comes standard with three shelves and a push handle. The top two shelves have three sides and an open front to slide in equipment with ease while also securing supplies. The bottom shelf has no sides for easy access to materials. We offer optional designs with additional guardrails and IV poles.

Our light duty carts have the same dimensions:

Width: 28 1/4" x Height: 35 3/4" x Depth: 16"



Top Shelf: Three Sided Shelf with Open Front Edge Middle Shelf: Three Sided Shelf with Open Front Edge Bottom Shelf: Open Shelf



Top Shelf: Three Sided Shelf with Open Front Edge Middle Shelf: Three Sided Shelf with Front Guardrail Bottom Shelf: Shelf with Four Guardrails



Top Shelf: Three Sided Shelf with Front Guardrail Middle Shelf: Three Sided Shelf with Front Guardrail Bottom Shelf: Shelf with Four Guardrails

Name	Light Duty Utility Cart	Light Duty Utility Cart with 5 Guardrails	Light Duty Utility Cart with 6 Guardrails
Model #	7534SS	7534SS-GR	7534SS-GR4
Item #	2427534002	2427534000	2427534004
Shipping Weight	30 lbs.	30 lbs.	30 lbs.

Options: IV pole attachments.

Custom Designs Available Upon Request

MEDIUM DUTY CARTS

Description: Our medium duty utility carts each come with three solid stainless steel shelves and travel on 4" swivel casters. Choose between our design that includes four sided shelves and width-wise handles or our design that includes open shelving and

depth-wise handles.



Dimensions: $27 \frac{1}{2}$ " x 40" x 18 $\frac{1}{8}$ "



Dimensions: 25" x 16 $\frac{1}{2}$ " x 35 $\frac{3}{4}$ "

Name	Side Rails Medium Duty Cart	End Rails Medium Duty Cart
Model #	7535SS	7536SS
Item #	2427535000	2427536000
Shipping Weight	47 lbs.	40 lbs.

Options: IV pole attachments.



Dimensions: 41 3/4" x 40" x 19 1/4"

HEAVY DUTY CART

Description: Our heavy duty utility cart comes with two solid stainless steel shelves with four sides each and travels on 5" swivel casters. Easy maneuverability with angled pushhandles.

Name	Heavy Duty Cart	
Model #	7537SS	
Item #	2427537000	
Shipping Weight	43 lbs.	

Options: IV pole attachments.

Custom Designs Available Upon Request



Accessories

IV POLE ATTACHMENTS

Description: Our IV pole attachment can fit any Blickman utility cart with ease.



Name	7534SS with IV pole attachment	7535SS with IV pole attachment	7536SS with IV pole attachment	7537SS with IV pole attachment
Model #	7534SS-IV	7535SS-IV	7536SS-IV	7537SS-IV
Item #	2427534001	2427535001	2427536001	2427537001
Shipping Weight	35 lbs.	52 lbs.	45 lbs.	48 lbs.

GUARDRAILS & SIDING

Contact Blickman Customer Service about adding additional siding or guardrails.







Open Shelf

Lipped Shelf

Shelf with 4 Guardrails

3 Sided Shelf

Custom Designs Available Upon Request



MR Conditional

Utility Carts

UTILITY CARTS

Description: Our MR conditional utility carts are made out of industry standard, fully welded, stainless steel construction. They feature two to three easy-access storage shelves and two ergonomic push handles. Optional IV pole attachments available for every model.



Name	Medium Duty End Rails Utility Cart	Medium Duty Side Rails Utility Cart	Heavy Duty Utility Cart
		MR Conditional Standard Models	
Item #	7536SS-MR	7535SS-MR	7537SS-MR
Model #	2497536000	2497535000	2497537000
Shipping Weight	40 lbs.	47 lbs.	43 lbs.
	MR C	Conditional Models with IV Attachments	
Item #	7536SS-IV-MR	7535SS-IV-MR	7537SS-IV-MR
Model #	2497536001	2497535001	2497537001
Shipping Weight	45 lbs.	52 lbs.	48 lbs.

Construction Specifications

Universal Dimensions	Light Duty: Width: 28 ¹ / ₄ " x Height: 35 ³ / ₄ " x Depth: 16"
Materials	Stainless steel Rubber bumper handles on 7534SS Carts
Gauge	Light Duty: 20GA & 18GA Medium/Heavy Duty: 16GA
Welding	All welding by spot and heliarc methods and exposed welds are ground smooth
Finish	Exposed surfaces are free from pit marks, weld seams, scale and are polished with a #4 polish to a satin finish
Casters	4" Swivel Casters

Shipping Information

All shipping weights are estimates and may be changed without prior notice. Specifications may be improved and dimensions changed without prior notice. Please contact customer service for estimated shipping costs at info@blickman.com.

Custom Designs Available Upon Request





Keeping patients warm and comfortable

Includes:

- Freestanding •
- **Under Counter** Recessed Accessories
- Tabletop



THE BLICKMAN **DIFFERENCE**

- We maintain the broadest line of warming cabinets in the healthcare market. Our blanket and fluid warming cabinets include freestanding, single chamber, dual chamber, tabletop, under-counter, and recessed models (including pass-through versions with doors on front and back).
- All of our models are offered with the choice of glass or stainless steel doors at the same price. Our glass
 door versions have wide windows designed to allow medical professionals to see all solutions and blankets
 at a glance, unlike the narrower porthole openings in competing models, so you never come up short of
 supplies in critical moments.
- Our warming cabinets allow for heat variances between 90° and 160° Fahrenheit. Temperatures remain within +/- 1° Fahrenheit of set point with virtually no leakage into adjoining champers or outside area.
- The warming fan is located at the top interior of the cabinet allowing for improved heat circulation (no hot/cold spots) and frees up the bottom shelf surface to store more blankets and/or solutions.
- Our touch screen and digital control panels allow O.R. managers to lock the warming cabinet settings with a key or passcode to ensure recommended temperature levels are maintained and that the facility meets all requirements.
- Unlike with cumbersome, "meat locker" handles on competing models, our warming cabinets allow for hands free operation. Our warming cabinets comes standard with self-closing, adjustable tension, magnetic gasket doors. If your arms are full, simply open the door with your elbow and then the door will slowly close on its own.

Key Product Features

Sloped top available for freestanding models to prevent dust collection

Dual chamber models enable you to heat solutions and blankets simultaneously in their respective recommended temperatures

Wide glass window designed to let medical professionals see all solutions and blankets at a glance

Free standing and recessed options available for tall models



Touch Screen Control Panel



Digital Screen Control Panel



Universal Controller Features:

- Visual/audible alarm
- Digital adjustments
- Calibration due date reminder
- · Graphical temp. trend view
- Set temp. "lock out"
- · Lost passcode view
- Temp. data recorder w/ real time clock

Self-closing, magnetic gasket, adjustable tension doors afford convenient, hands free opening and closing

Freestanding

Description: Our freestanding cabinets come in either single or dual chamber configurations with open adjustable shelf design for flexibility and expanded use. The doors come in either glass or stainless steel and lock to secure contents. Each unit comes standard with three perforated stainless steel shelves and integrated casters for ease of mobility.

Dimensions

Width: 30" x Height: 74 $\frac{1}{2}$ " x Depth: 26 $\frac{5}{8}$ "









Electrical: Amps: 7.5 | Total Wattage: 750

Electrical: Amps: 15 | Total Wattage: 1500

Name Single Chamber			Dual C	hamber	
			Glass Doors		
Controller Type	Touch Screen	Digital Screen		Touch Screen	Digital Screen
Model #	8921TG	7921TG		8924TG	7924TG
Item # Shipping	14B8921243	14B7921243		14B8924243	14B7924243
Weight	437 lbs.	437 lbs.		500 lbs.	500 lbs.
			Solid Doors		
Controller Type	Touch Screen	Digital Screen		Touch Screen	Digital Screen
Model #	8921TS	792ITS		8924TS	7924TS
Item # Shipping	14B8921200	14B7921200		14B8924200	14B7924200
Weight	437 lbs.	437 lbs.		500 lbs.	500 lbs.

Options: No Casters, sloped tops

Recessed

Description: Our recessed warming cabinets are designed for built-in installation. Single chamber glass door, solid door, and passthrough models available. Three perforated stainless steel shelves come standard in these models.

Dimensions

Width: 30" x Height: 60" x Depth: 20 5/8" *Electrical*: Amps: 7.5 | Total Wattage: 750

Glass Door



Solid Door



Passthrough



Name	Glass	Door	Solid	Door	Passtl	hrough
Controller Type	Touch Screen	Digital Screen	Touch Screen	Digital Screen	Touch Screen	Digital Screen
Model #	SW40TG	SW30TG	SW40TS	SW30TS	SW40TPTGS	SW30TPTGS
Item # Shipping	14BSW40243	14BSW30243	14BSW40200	14BSW30200	14BSW40278	14BSW30278
Weight	320 lbs.	320 lbs.	320 lbs.	320 lbs.	245 lbs.	245 lbs.

Options: Dual chambers, additional shelves.

Note: Passthrough models are available in the following configurations: glass/glass, solid/solid, glass/solid, or solid/glass.

Under Counter

Description: Our under counter warming cabinet is suited for use in limited spaces. It comes with either glass or stainless steel doors, and comes standard with two adjustable perforated stainless steel shelves.

Dimensions

Width: 30" x Height: 35 1/2" x Depth: 26 5/8"

Electrical: Amps: 7.5 | Total Wattage: 750





Name Single Chamber

Glass Door		Solid	Door	
Controller Type	Touch Screen	Digital Screen	Touch Screen	Digital Screen
Model #	8922TG	7922TG	8922TS	7922TS
Item # Shipping	14B8922243	14B7922243	14B8922200	14B7922200
Weight	260 lbs.	260 lbs.	260 lbs.	260 lbs.

Options: Sloped tops

Tabletop

Description: Our tabletop warming cabinets are available in a standard and compact size. They can sit on an existing countertop or on a Blickman mobile stand specifically designed to support our tabletop warmers. They come with either glass or stainless steel doors. One adjustable perforated stainless steel shelf is included in the standard model only.

Dimensions: 30" x x 24 1/2" x 26 5/8" **Electrical:** Amps: 7.5 | Total Wattage: 750





Dimensions: 24" x 24 1/2" x 20 5/8" **Electrical:** Amps: 3.5 | Total Wattage: 350





Name	Stand	dard		Con	npact
			Glass Doors		
Controller Type	Touch Screen	Digital Screen		Touch Screen	Digital Screen
Model #	8925TG	7925TG		8927TG	7927TG
Item # Shipping	14B8925243	14B7925243		14B8927243	14B7927243
Weight	215 lbs.	215 lbs.		161 lbs.	161 lbs.
			Solid Doors		
Controller Type	Touch Screen	Digital Screen		Touch Screen	Digital Screen
Model #	8925TS	7925TS		8927TS	7927TS
Item # Shipping	14B8925200	14B7925200		14B8927200	14B7927200
Weight	215 lbs.	215 lbs.		161 lbs.	161 lbs.





Description: Our warming cabinet mobile stand allows for our tabletop warming cabinets to become mobilized for transport around your medical facility. These tables are designed to fit both our standard or compact tabletop warming cabinet dimensions without taking up any extra floor space. It includes three sides to secure warming cabinet placement and an optional shelf with four guardrails is available for additional storage.

For transportation, wrap warming cabinet cable around cleat on the back of the mobile stand.

Table Dimensions

Width: 30" x Height: 34" x Depth: 26"



Cleat in the back to coil power cord during transport.





Name	H-Brace	Shelf
Model #	3026SSH	3026SS
Item #	0123026100	0113026000
Shipping Weight	40 lbs.	50 lbs.

Accessories & Power Consumption

ACCESSORIES



Item	Description	Shipping Weight
14B3018000	Sloped top (30" x 8" x 18") for warming cabinet SW40	28 bs.
14B3026000	Sloped top (29 ⁷ / ₈ " x 10" x 26") for warming cabinets: 8921, 8922, 8924, 8925	38 lbs.



Item	Description	Shipping Weight	
14B8924SHF	Warming Cabinet Shelf Kit 30" x 26" with 1" Diameter holes	78 lbs.	
	for 8921, 8922, 8924, 8925		
14B8927SHF	Warming Cabinet Shelf Kit 20" x 17 1/2" for 8927	28 bs.	
14BSW40SHF	Warming Cabinet Shelf Kit 30" x 18" for SW40	38 lbs.	

POWER CONSUMPTION

Model#	Description	Cubic Feet	KW/Hr.	BTU/Hr.
8921	Single Chamber Tall Warming Cabinet	22.02	0.55	1878
8922	Under Counter Warming Cabinet	9.02	0.35	1195
8924	Dual Chamber Tall Warming Cabinet	18.58	0.78	2664
8925	Standard Tabletop Warming Cabinet	5.05	0.29	990
8927	Compact Tabletop Warming Cabinet	2.87	0.38	1298
SW40	Recessed Warming Cabinet	13.22	0.47	1605



Construction Specifications

	Construction Specifications
Standard Dimensions	Free Standing: 30" x 74 1/2" x 26 5/8" Recessed: 30" x 60" x 20 5/8" Table Top: Standard: 30" x 26 5/8" x 24 1/2"; Compact: 24" x 20 5/8" x 24 1/2" Under Counter: 30" x 60" x 20 5/8"
Materials/ Gauge	304 stainless steel 1/4" clear double panel, insulated, tempered sheet glass Perforated Stainless Steel Shelves
Welding	All welding by heliarc method and exposed welds are ground smooth
Finish	Exposed surfaces are free from pit marks, weld seams, scale and are polished with #4 polish to a satin finish
Electrical	Available in standard 100/120 VAC input or 220/240 VAC Units are FDA registered and conform to standards for Safety Medical Electrical Equipment UL 606011 (US) and CSA C22.2 #601.1M90 (Canada) Units are cord connected with a hospital grade cap
Casters	Tall Warming Cabinets: Hidden 3" ball bearing polyolefin casters with zinc finish

Notes •

- TOUCHSCREEN CONTROLLER: LCD display shows actual temperature and set point for each chamber at all times. Control panel components include an On/Off Switch, Keyed Control Temperature Lockout USB port and Touchscreen User Interface.
- DIGITAL CONTROLLER: Dual display shows actual temperature and set point for each chamber at all times (each chamber in the double door warmer includes its own digital control system). Control panel components include an ON/OFF Button, UP/DOWN Temperature Touch pad, Alarm Light, Heater Light Indicator, Set/Mute Touch
- TEMPERATURE: 90° to 160° Fahrenheit. There is a heater light indicator, digital temperature indicators and an over temperature alarm indicator if the actual temperature exceeds set temperature by 10° Fahrenheit. All units have a visual and audible alarm as well as an automatic safety switch. The units are furnished with tubular heater assemblies and a circulating fan for even heat distribution. From a cold start, a fully loaded cabinet will be heated to the selected temperature in four to six hours.

Shipping Information

All shipping weights are estimates and may be changed without prior notice. Specifications may be improved and dimensions changed without prior notice. Please contact customer service for estimated shipping costs at info@blickman.com.

Contact Us

To place an order or request a quote, go to our website:

https://www.blickman.com/order-form https://www.blickman.com/request-a-quote-form

If you have any questions or would like to speak with our customer service department:

Toll Free: (800) 247-5070 **Phone:** (973) 330-0557

Email: info@blickman.com **Website:** www.blickman.com

If you require tech support, please visit our website at www.blickman.com/tech-support or email us directly at techsupport@blickman.com.

OUR PROMISE

For over 125 years, Blickman has been known throughout the healthcare industry for exceptional product quality. We pride ourselves on providing high-end and personalized service to ensure 100% customer satisfaction. If you are not completely satisfied for any reason, we will make it right. All of our equipment is backed by our Limited Lifetime Warranty, so you can be certain that you have made the very best choice for your healthcare facility.



www.blickman.com









Corporate Office & Headquarters

1099WallStreetWest Suite250 Lyndhurst,NJ 07071