

BLICKMAN



Solution Stands/ Receptacles

Roll, collect, dispose

Includes:

- Circulator Station
- Kickbucket
- Receptacle
- Solution Stands



THE BLICKMAN DIFFERENCE

- Our high quality, heavy-gauge, noise-reducing, solution stands and receptacles ensure your solutions and waste products are easily collected and disposed of in these easy-to-move, easy-to-clean product.
- Our solution stands and receptacles move effortlessly and quietly on 2" swivel casters, and are made out of our high-quality, easy-to-clean 304 stainless steel with MR conditional options available for both.
- Our solution stands are 33" high for compatibility with surgical table heights and are in the safe zone for reduced back stress, and our receptacles are low to the floor to pass under surgical tables and make it easy to drop waste into without reaching around the table.
- Our basins and pails are removable for both the solution stands and receptacles that are held in place by rubber bumpers that also help prevent noise from metal on metal contact.

Key Product Features



Circulator Station

Promote safety, communication, and efficiency in the OR.

Description: This central station includes a communication board, a utility shelf, and sponge count bags. The communication board improves patient-centered care, as it can be used for displaying critical medical information to the entire surgical team, such as specific information about the patient, the surgery safety checklist (SSC), doctor's requests, and more. This improved communication results in decreased sentinel events and human factor errors. The multi-utility surface acts as a rotational shelf to hold intra-operative medications, small counted items, and specimens that are visible to all surgical team members and hangs up to five standard sponge count bags.



20 7/8"W x 18 1/4"L x 71-91"H

CIRCULATOR STATION

Name	Model #	Item #	Ship Weight
Circulator Station	2686SS	2732686000	41 lbs.

Receptacles

Description: Our stainless steel kickbucket consists of a continuous encircling rubber bumper, which helps prevent tipping, wall scuff marks, and noise, a stainless steel cross brace, which provides excellent protection and a secure hold of a heavy gauge 13 qt. pail, and roller blade casters for mobility.

Parts Also Sold Separately



Kickbucket

13 qt. Pail



#KB-1



#7766NP

Parts Also Sold Separately



Receptacle

8.5 qt. Solution Basin





#SB-1



#8766NP

KICKBUCKET

Name	Model #	Item #	Ship Weight
Lenox Kickbucket with Pail 	7766SS	0817766000	12 lbs.
Lenox Kickbucket with Pail 	7766MR	0817766066	12 lbs.

ORDER SEPARATELY

Lenox Kickbucket Frame Only	7766NP	0817766001	9 lbs.
13 qt. Pail Only	KB-1	9701920000	4 lbs.

RECEPTACLE

Name	Model #	Item #	Ship Weight
Sponge Receptacle with Basin	8766SS	0828766000	11 lbs.

ORDER SEPARATELY

Sponge Receptacle Frame Only	8766NP	0828766001	9 lbs.
8.5 qt. Solution Basin Only	SB-1	9737370120	3 lbs.

Solution Stands

Description: We offer two types of stainless steel solution stands: a single basin model and a double basin model. Our solution stands are 33" tall offering compatibility with surgical table heights and are in the safe zone for reduced back stress. Choose between a lower shelf or an H-brace for added storage and support. MR conditional models available.









Dimensions:

Width: 15" x Height: 33" x Depth: 15"



Dimensions:

Width: 29" x Height: 33" x Depth: 15"

Name	Single Basin Solution Stand		Double Basin Solution Stand		
			Shelf		
					
	Model #	7807SS	7807MR	7808SS	7808MR
	Item #	0717807000	0717807066	0727808000	0727808066
	Shipping Weight	20 lbs.	20 lbs.	33 lbs.	33 lbs.
		H-Brace			
	Model #	7807SS-HB	7807MR-HB	7808SS-HB	7808MR-HB
	Item #	0717807100	0717807166	0727808100	0727808166
	Shipping Weight	19 lbs.	19 lbs.	24 lbs.	24 lbs.

Solution Basin Also Sold Separately



Name	Model #	Item #	Ship Weight
8.5 qt. Solution Basin Only	SB-1	9737370120	3 lbs.



Solution Stands

Construction Specifications

Construction Specifications

Universal Dimensions	Solution Stands: 15", 29" x 33" x 15" Solution Basin: 8.5 qts. Pail: 13 qts.
Materials	304 stainless steel
Gauge	Solution Stand Frame: 1" diameter tube, 16GA Solution Stand Shelf: 20GA Pail: 22GA; Basin: 14GA
Welding	All welding by heliarc method and exposed welds are ground smooth
Finish	Exposed surfaces are free from pit marks, weld seams, scale and are polished with #4 polish to a satin finish
Casters	2" swivel casters
Notes	<ul style="list-style-type: none"> MR conditional options available where indicated

Shipping Information

All shipping weights are estimates and may be changed without prior notice. Specifications may be improved and dimensions changed without prior notice. Please contact customer service for estimated shipping costs at info@blickman.com.

Toll Free: (800) 247-5070
Local: (973) 330-0557



www.blickman.com
info@blickman.com

1099 Wall Street West, Suite 250A, Lyndhurst, NJ 07071