

Foot Stools

Safe and sturdy to step up or down

Includes:

- 8" High Foot Stools
- 9" High Foot Stools
- Two Step Foot Stools
- EZ Stacking Stools + Dolly



THE BLICKMAN DIFFERENCE

- Whether chrome or stainless steel, our foot stools are high quality, durable, and safe.
- Every foot stool includes a non-slip rubber tread for secure footing and non-slip rubber tips on the stools' feet to secure floor placement.
- For heavy loads, some models include H-brace reinforcement.
- For additional user stability, some models include handrails or allow for the addition of a handrail.
- MR conditional options are also available.

Key Product Features





One Step

8" HIGH FOOT STOOLS

Description: Our MR conditional approved Kent foot stools include one step 8" off the ground. They are industry standard, fully welded, stainless steel construction with non-slip rubber treads and feet tips. 18" and 30" width sizes are available. Depth is 12". Include a 41" tall handrail for additional support.

8" stools have the same dimensions:

Width: 18" or 30" x Height: 8" x Depth: 12"



Name	18"W	30″W	
	With Handrail		
Model #	7757MR-HR	7761MR-HR	
Item #	1017757100	1017761100	
Shipping Weight	20 lbs.*	30 lbs.*	
	With	out Handrail	
Model #	7757MR	7761MR	
Item #	1017757000	1017761000	
Shipping Weight	11 lbs.	20 lbs.	



One Step

9" HIGH FOOT STOOLS

Description: Our 9" high foot stools come in either chrome or stainless steel and consist of one step 9" off the ground with non-slip rubber treads and feet tips. A 35 3/4" tall handrail for additional support is available for the chrome stool.

All 9" stools have the same dimensions:

Width: 14 1/4" x Height: 9" x Depth: 11 3/8"





Name	Stainless Steel Stool	Chrome Stool	
	With Handrail		
Model #	1261	1251	
Item #	1011261000	1011251000	
Shipping Weight	12 lbs.*	11 lbs.*	
	Without	t Handrail	
Model #	1260SS	1250	
Item #	1011260000	1021251640	
Shipping Weight	8 lbs.	7 lbs.	



Foot Stools

Two Step

TWO STEP FOOT STOOLS

Description: Our MR conditional approved two step foot stools are industry standard, fully welded, stainless steel construction with H-braces for additional reinforcement, non-slip rubber treads for secure footing, and feet tips for firm floor placement.



Step Surface Areas Top Step: 36" x 12"

Top Step: 36" x 12" Bottom Step: 36" x 24"

Heights of Steps

Top Step: 15 3/4" Bottom Step: 8"

Name	2-Step Foot Stool	
	With Handrails (41" high)	
Model #	7764MR-HR	
Item #	1012041096	
Shipping Weight	65 lbs.	

Two Step Foot Stool with Side Handrails



The Donnelly with Handrail

Step Surface Areas

Top Step: 24" x 8" Bottom Step: 24" x 16"

Heights of Steps

Top Step: 15 3/4" Bottom Step: 8"

Name	The Donnelly	
	With Handrails (41" high)	
Model #	7763MR-HR	
Item #	1017763100	
Shipping Weight	40 lbs.*	
	Without Handrail	
Model #	7763MR	
Item #	1017763000	
Shipping Weight	27 lbs.	



EZ Stacking

OUR EZ Stacking Stool allows you to:

- ✓ Create wider platforms by locking side-to-side
- ✓ Create higher platforms by stacking top-to-bottom
- ✓ Store efficiently and transport easily with our corresponding dolly

EZ STACKING STOOLS



Description: Our EZ stacking stools and dolly are made with our durable stainless steel. The stools are covered with a non-slip rubber tread and feet tips.



EZ Stool Dimensions: 18 ⁷/₈" x 6 ¹/₄" x 13 ⁷/₈"



Dolly Dimensions: $21 \frac{1}{2}$ " x 39" x $16 \frac{1}{2}$ "

Name	EZ Stacking Stool	Dolly
Model #	8862SS	8862SS-D
Item #	1018862000	1018862300
Ship Weight	14 lbs.	20 lbs.

Construction Specifications

Standard 1 step: 6 1/4" - 9" high

Dimensions 2 step: Bottom step: 8" high, Top step: 15 3/4" high

Materials Chrome

304 stainless steel Black non-slip rubber

Gauge EZ Stacking Stools: 16GA platform, 14GA legs

All Stainless Steel Stools: 16GA (Exception: 7757MR-HR and 7761MR-HR

have 18GA handrails)

Welding All welding by heliarc method and exposed welds are ground smooth

Finish Exposed surfaces are free from pit marks, weld seams, scale and are

polished with #4 polish to a satin finish

Casters Dolly: 3" swivel casters with no brakes

Notes • MR conditional foot stools are available where indicated

Note for:

Parts shipped separately to avoid damage and lower costs. Shipping weight with asterisk is the combined shipping weight of separate parts. Some assembly required.

1251 Stool

Shipping Information

All shipping weights are estimates and may be changed without prior notice. Specifications may be improved and dimensions changed without prior notice. Please contact customer service for estimated shipping costs at info@blickman.com.

Toll Free: (800) 247-5070 Local: (973) 330-0557



www.blickman.cominfo@blickman.com